
# TheRoyal

DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at [denroyal.com/symbols](http://denroyal.com/symbols) or by contacting customer service.

## EN: FOLEY CATHETER WITH TEMPERATURE SENSOR

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## FOLEY CATHETER WITH TEMPERATURE SENSOR

MR	MR CONDITIONAL
	<b>MEDICAL DEVICE</b>
	<b>DO NOT RESTERILIZE</b>
	<b>DO NOT REUSE</b>
	<b>STERILIZED WITH ETHYLENE OXIDE</b>
	<b>NOT MADE WITH NATURAL RUBBER LATEX</b>
<b>RX ONLY</b>	<b>FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.</b>

**IMPORTANT INFORMATION**  
Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

**INTENDED USE**  
The DeRoyal Foley Catheter with Temperature Sensor is to be used for drainage of the urinary bladder and simultaneous monitoring of core body temperature.

**INDICATIONS**  
Foley catheterization is recommended when continuous bladder drainage is required in the absence of voluntary micturition, after surgery or trauma of the pelvic organs, in the presence of obstruction or paralysis, and as a means to evaluate urine output in patients having fluid management problems.

**DESCRIPTION**  
The Foley Catheter with Temperature Sensor consists of a temperature-sensing wire set secured inside a silicone catheter. The catheter's balloon inflation system is terminated with a luer tip syringe-actuated, self-sealing valve. The device is available in 8 - 18 French size. The French size and balloon volume are designated on the catheter.

The wire set is secured within the lumen of the foley catheter. The temperature sensor is located near the distal tip of the catheter. The connector at the proximal end interfaces with DeRoyal approved cables for connection to a monitoring unit. The temperature monitoring sensor is available in 400 or 700 series configurations for connection to a compatible monitoring unit. Consult the user manual for the patient monitor to ensure compatibility.

**WARNING**  
The patient should be routinely monitored for urine output to minimize risk of urine retention.  
Use aseptic technique in patient preparation and catheter care.  
Follow accepted medical techniques during removal of catheter to minimize risk for patient infection or injury.  
Carefully inspect the packaging and device for damage or defect. **DO NOT** use if the packaging is open or the device is damaged.

**DO NOT** use device on a patient with weak bladder walls.  
Chronic catheter use increases the risk of bladder retention, urinary tract infection or other adverse event. Minimize catheterization duration following accepted clinical practice. If resistance to flow in the catheter develops, then the catheter should be changed.  
Dispose of device according to facility protocol.  
**DO NOT** store device in extreme temperature or a moist environment.  
This device is for single patient use only. It is not intended for disinfection and subsequent re-use, which may result in device failure and/or create the risk of contamination.  
**DO NOT** use an expired device.  
Prior to removing the catheter, deflate the balloon by inserting the luer tip of an empty syringe into the valve and aspirating the solution completely. If the balloon does not collapse after attempting this step, then sever the

valve arm. If balloon still does not deflate, call physician.

### ADVERSE REACTIONS

The following have been reported to be associated with the use of foley catheter and/or temperature sensors, especially during prolonged use:

- bladder spasms
- electric burns due to aberrant electrosurgical current
- encrustations with urinary salts
- kidney or bladder stones
- septicemia
- urethral discharge
- urethritis
- urethrocutaneous fistulae
- urinary tract infection or irritation

### CAUTIONS

- Take proper care to not damage the insulation prior to or during use.
- Prior to insertion, test the balloon inflation and deflation. The balloon should be fully deflated before insertion.
- Follow the recommended balloon inflation volume printed on the product and use sterile water for inflation.
- Select the appropriate size catheter using expert clinical and medical judgment.
- Lubricate catheter with a suitable water soluble lubricant prior to insertion
- Deflate the balloon prior to removal of the catheter.
- DO NOT** use a needle to puncture the catheter wall for aspiration. Leaks may occur.
- DO NOT** intertwine the cables, especially the monitor cables with the electrosurgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electrosurgical activations. Unusual temperature readings should be checked.
- The catheter is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

### DIRECTIONS FOR USE

- Use aseptic technique in patient preparation and catheter care.
- Verify the compatibility of the catheter, interface cable, and patient monitor before use.
- Remove the Foley Catheter with Temperature Sensor from its sterile package.
- Lubricate the catheter with a suitable water soluble lubricant and insert using acceptable medical techniques.
- Make sure that the catheter is draining properly before inserting the balloon.
- To inflate the balloon, insert a luer tip syringe filled with the recommended amount of sterile water into the valve. Depress the plunger. The valve will seal when the syringe tip is withdrawn.
- Attach drainage tubing and bag.
- Secure the catheter, drainage tubing and bag in place by following accepted medical technique.
- Connect the temperature probe to the appropriate reusable cable from DeRoyal. Connect the cable to the patient monitor. Secure the cable using the drape clip.
- Follow the directions for use of the patient temperature monitor instrument.
- Prior to removing the catheter, deflate the balloon by inserting the luer tip of an empty syringe into the valve and aspirating the solution completely.
- If balloon still does not collapse after attempting to perform #10 above, sever the valve arm. If balloon still does not deflate, call physician.
- Remove the catheter following accepted medical techniques.
- Disconnect the probe (at the connector) and discard probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.
- Dispose of probe according to facility protocol.

### MR COMPATIBILITY

The MR Catheter with Temperature Sensor was determined to be MR Conditional according to the terminology specified in the American Society

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For Testing and Materials (ASTM) International designation: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment, and other related standards.

This Foley Catheter with Temperature Sensor is safe for use with MR imaging, when utilized in compliance with the instructions located here: [www.denroyal.com/foley-catheter-mr-compatibility](http://www.denroyal.com/foley-catheter-mr-compatibility) or by calling: 1-800-251-9864. **WARNING: USE OUTSIDE OF THESE PARAMETERS COULD RESULT IN SERIOUS INJURY OR DEATH.**

### OPERATING SPECIFICATIONS

- Rated output range: 25° to 45° Celsius
- The reference body site is the core body temperature.
- The measuring site is the Bladder.
- Accurate to +/- 0.2° Celsius
- Time response:
  - The temperature probe has a direct mode of operation.
  - The maximum heating transient time is 1 minute
  - The maximum cooling transient time is 30 seconds

### STORAGE AND TRANSPORT CONDITIONS

	<b>KEEP DRY</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>
	<b>DO NOT USE IF PACKAGE IS DAMAGED</b>
	<b>RELATIVE HUMIDITY: 85% NON-CONDENSING</b>
	<b>STORAGE TEMPERATURE: 25° C TO + 55° C</b>
	<b>OPERATING TEMPERATURE: 25° C TO 45° C</b>

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

### WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES. INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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