

Comfort Cool™ Thumb CMC Restriction Splint

Single Patient Use Only

Rx Only

IMPORTANT INFORMATION

Please read all instructions, warnings, and precautions before use.

INTENDED USE



To provide thumb CMC joint support and light compression for CMC joint arthritis and/or instability, or for discomfort associated with repetitive motion. May be helpful with tenosynovitis of the thumb extensors/abductors.

PRECAUTIONS

Consult your clinician immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this product.

INSTRUCTIONS FOR USE

1. Slide thumb into the thumb sleeve. Wrap the splint around the wrist and secure the hook and loop fastener comfortably on the dorsum (top) of the hand.
2. Pull the CMC strap firmly from the base of the thumb, attaching it to the curved hook fastener on the thenar (palm) eminence.



3. Gently pull the strap through the thumb web space and secure it to the hook fastener on the dorsum (top) of the splint. Trim the strap if needed.



4. To allow greater IP motion, cut curve along volar thumb edge. If the thumb MP joint support is not desired, cut off the thumb sleeve below the thumb web seam. Be sure to trim just above one of the stitched reinforcements. Trim next to the distal edge of the curved hook fastener on the thenar no closer than 1/4" (6.4mm).



SIZING

Made for hands with enlarged thumbs, plus sizes are designed one size wider in the thumb and a half size bigger in the hand.

MP CIRCUMFERENCE

Small	6" – 7"
Medium	7" – 8"
Large	8" – 9"

CARE INSTRUCTIONS

Hand wash with mild soap in cool water and air dry for best results. Secure hook and loop during washing.

Comfort Cool™ is a trademark of North Coast Medical.



Distributed by:
DeRoyal
200 DeBusk Lane, Powell, TN 37849 U.S.A.
888.938.7828
865.938.7828

©2004 DeRoyal. All Rights Reserved.
Reprint #1-348 Revision 01
Revised 2/04