SHOULDER ABDUCTION DEVICE

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SCHULTERABDUKTIONSORTHESE

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Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (601) 865.938.7828
www.derooyal.com

EMERGÉ EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands

Part#0-038G
Revised 2/2021

©2021 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.

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SHOULDER ABDUCTION DEVICE

	SINGLE PATIENT USE
	NON-STERILE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury or re-injury inherent with the use of any brace.

INTENDED USE

The DeRoyal® Shoulder Abduction Device is intended for patients recuperating from rotator cuff repair, Bankart repairs, soft tissue repairs, shoulder sprains and strains, acute shoulder instability and reconstruction.

CAUTIONS

- A physician or properly licensed practitioner (a “prescriber”) who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber. Use only as directed.
- Prolonged usage can result in muscle atrophy in the user. The user may require (1) a wearing period from the brace to build muscle strength (2) regular physical therapy exercises to maintain the strength of the shoulder (3) a limitation on duration of use or (4) some combination of these.
- Inspect brace for damaged or missing components before use.
- Discontinue use and consult your prescriber if the brace or its components break or become damaged.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.
- Take special care if the user is diabetic or has poor circulation. These users may have decreased skin sensitivity and are at greater risk for poor peripheral circulation and pressure injuries.

INSTRUCTIONS FOR USE

- Place cushion under affected arm as high as possible.
- Attach 4” (10 cm) wide hook and loop backed panel of chest strap to side of cushion with top

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end of panel slightly above lower angle of cushion. **FOR RIGHT ABDUCTION:** Place short end of strap toward the back and close with hook and loop. Wrap long end of chest strap across chest, under arm and across the back. **FOR LEFT ABDUCTION:** Place short end of strap toward front and close with hook and loop. Wrap chest strap across the back, under arm and across chest. Attach closure to part of the chest strap already in place.

- In both applications, pull chest strap tightly before attaching closures, keeping the cushion as high as possible. If cushion can be pushed downward with hand pressure, then the strap is not tight enough.
- Wrap upper arm cuff around arm above elbow and attach to cushion to help hold arm in place.
- Wrap wrist cuff around wrist and upper part of hand and fasten with hook and loop. Be sure the strap hangs down from the bottom of the cuff. Position hook and loop backed panel of cuff onto cushion. Bring wrist strap up and across unaffected shoulder and attach hook and loop panel to chest strap in back. If necessary, adjust the buckle and cut off the excess strap.

CARE INSTRUCTIONS (IF APPLICABLE)

Hand wash with mild detergent in lukewarm water. Air dry prior to wearing.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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