

anchoring them at the hook strips found on the back of the lateral panels.


- Attach the gray loop material on the left side of the brace belt to the gray hook material on the anterior panel. Ensure the oval shape of the belt flap is evenly aligned to the oval shape of the hook on the anterior panel. **NOTE:** Ensure both Boa® Fit System dials are in the released (pulled-out) position.
- Close the belt by attaching the right side of the belt evenly over the left side of the belt and anterior panel. Ensure the oval shape of the right side of the belt is evenly aligned to the oval shape of the left side of the belt and the lateral panels.
- Retrieve the Sternal Pad straps from behind the posterior panel and attach the straps via the quick release buckles (2) on the sternal pad plate. **NOTE:** The straps can be positioned either over the shoulder or under the arm for desired comfort.
- Tighten the brace belt by pushing in on the Boa® Fit System dials on each side of the brace and turning evenly to apply compression. Each dial should be turned in the same direction (clockwise) and at the same time to achieve even compression. Tighten until desired compression is achieved. **NOTE:** Ensure the orthosis is in a midline position with the spine and that it has not migrated upwards or downwards.
- Tighten the Sternal Pad straps at the quick release buckles (2). **NOTE:** Straps may be trimmed, if needed.

## PANEL MODIFICATION FOR BRACE CUSTOMIZATION

**NOTE:** See customization worksheet for further details and documentation of brace customization for patient.

- The posterior panel of the brace can be modified using a regulated heat gun. Remove the Posterior Thoracic Extension (PTE) from the top of the posterior panel (see below). Position the posterior panel in the appropriate lordotic curvature, apply heat, and allow to cool in that position. Modification is only recommended above the B and D ports as indicated on the inside of the posterior panel. The adhered hook may be removed during modification and reapplied once material has cooled. Let panel cool completely before applying brace to patient. **NOTE: DO NOT** attempt to heat mold around the PTE engagement area as this may compromise the PTE engagement mechanism.
- The lateral panels on each side of the back panel can be trimmed or removed for a more customized and comfortable fit.

## POSTERIOR THORACIC EXTENSION (PTE) REMOVAL

If the patient progresses or the PTE is not needed, it can be removed by lifting on the middle tab at the bottom of the PTE and pulling upwards. Ensure that the two feet of the PTE are disengaged at the same time as the lift tab to completely remove PTE from the brace.

## CARE INSTRUCTIONS

Remove panels from the brace and clean all soft goods by hand washing with mild soap and let air dry completely prior to reapplication. Anterior and posterior panels can be cleaned by wiping down with a disinfectant wipe.

## STORAGE AND TRANSPORT CONDITIONS

	<b>KEEP DRY</b>
	<b>KEEP AWAY FROM SUNLIGHT</b>

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

## WARRANTY

DEROYAL® PRODUCTS ARE WARRANTED FOR NINETY (90) DAYS FROM THE DATE OF SHIPMENT FROM DEROYAL AS TO PRODUCT QUALITY AND WORKMANSHIP. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

Manufactured by:  
DeRoyal Industries, Inc.  
200 DeBusk Lane  
Powell, TN 37849 USA  
888.938.7828 or (001) 865.938.7828  
www.deroyal.com

PART#0-1990 | REVISED 3/2023

©2023 DEROYAL INDUSTRIES, INC.

All Rights Reserved. MAXAlign, DeRoyal, and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

Boa is a registered trademark of Boa Technologies Inc.

## MAXALIGN® TLSO-S SPINE BRACE

	<b>SINGLE PATIENT USE</b>
	<b>NON-STERILE</b>
	<b>NOT MADE WITH NATURAL RUBBER LATEX</b>
<b>RX ONLY</b>	<b>FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.</b>

## IMPORTANT INFORMATION

Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of product.

## INTENDED USE

The DeRoyal® MAXAlign® TLSO-S Spine Brace powered by the Boa® Fit System is intended to provide support and immobilization of the spinal column in the sagittal, coronal and transverse planes. The device is also intended to aid in relief and recovery following spinal surgery or an injury.

## INDICATIONS

Post-operative thoracic/lumbar fusion, laminectomy, discectomy, compression fractures, degenerative disc disease, osteoporosis, single column instability, and facet syndrome.

## ⚠ CAUTIONS

- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.
- Inspect product for damage and/or missing components before use.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.**
- Prolonged usage can result in muscle atrophy in the user.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace.

## INSTRUCTIONS FOR USE

### SIZING

The DeRoyal MAXAlign® TLSO-S features 6 sizing zones:

**NOTE:** Brace sizing is determined by circumferential waist measurement taken at the navel.

	ZONE 1	ZONE 2	ZONE 3	ZONE 4	ZONE 5	ZONE 6
<b>MEASUREMENT</b>	24"–30"	30"–36"	36"–42"	42"–48"	48"–54"	54"–60"

**NOTE:** For best results and ease of adjustment, brace sizing should be conducted on a large bed or table.





## BRACE SIZING ADJUSTMENT:

- Prior to application, release the sternal pad straps from the sternal pad bar by pressing the quick release buckles.
- Pull out both Boa® Fit System dials to release lacing and




release the front hook and loop closure to open the brace belt completely and remove the sternal pad bar/anterior panel assembly from the brace and set aside for now.

- Place the back of the brace flat with the outside of the back panel facing upward. Ensure the left and right belt arms are spread outward evenly and are non-impinged.
- Unfasten the back-panel hook closures and flip over the back-panel and belt arms to reveal the belt sizing adjustment mechanism.
- Pull up on the two (2) belt security hook flaps located on each side of the adjustment mechanism until they reach a 45-degree angle upward from the belt.
- Place one hand under the Belt Pull Strap near the hand indicator and use the other hand to pull up on the Belt Pull Strap until the appropriate size zone is reached on the left and right side of the belt. **NOTE:** Ensure that the sizing zones are even on both sides. Reapply the belt security hook flaps to secure the belt.
- Align the Belt Pull Strap with the middle of the belt adjustment mechanism panel, ensure excess belt length is even and press down firmly. Fold the back-panel and belt over and reattach the back-flap hook closures at the top of the back panel.

## BRACE APPLICATION:

- To apply the brace in the supine position, roll the patient on their side. Tuck the side closure behind the posterior panel and slide the closure and lateral panel under the patient's lower back so the panels are accessible and can be retrieved on the opposite side.
- The inferior (lower) portion of the posterior panel should be positioned just above the patient's coccyx.
- Roll the patient back into the supine position.
- Next, retrieve the sternal pad bar/anterior panel assembly and adjust the sternal pad bar height by pinching the sides of the adjustment mechanism on the anterior panel and releasing when appropriate height is achieved. **NOTE:** Ensure the locking peg of the adjustment mechanism is engaged. The sternal pad bar is pre-bent to accommodate patient anatomy. Bar and plate can be bent further to adjust to patient anatomy as needed.
- Place the anterior panel in the appropriate position just above the pubic bone.
- Ensure the inside lateral panels are properly aligned and even with the brace belt by

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