



- retrieved after rolling the patient back into the supine position.
- The inferior (lower) portion of the posterior panel should be positioned just above the patient's coccyx.
 - Roll the patient back into the supine position.
 - Retrieve the brace belt closure from under the patient and leave extended out to the side.
 - Retrieve the DLE straps from above the patient's shoulders and anchor them with the hook and loop straps to the posterior panel in zone A, B, or C (labeled on the posterior panel) but do not tighten at this time. **NOTE:** For greater shoulder restriction, position A should be used.
 - Once the DLE straps are secured, continue by wrapping the left side of the brace belt around the body and position just above the pubic bone. **NOTE:** Ensure both Boa® Fit System dials are in the released (pulled-out) position.
 - Close the belt by attaching the right side of the belt evenly over the left side of the belt. Ensure the oval shape of the right side of the belt is evenly aligned to the oval shape of the left side of the belt.
 - Tighten the brace belt by pushing in on the Boa® Fit System dials on each side of the brace and turning evenly to apply compression. Each dial should be turned in the same direction (clockwise) and at the same time to achieve even compression. Tighten until desired compression is achieved.
 - Tighten the DLE straps using the excess strap under the patient's arms. Tighten to appropriate comfort and support. Excess floating strap may be tucked into the belt, if desired.
 - Position and secure the chest strap as needed along the tubing via the slide buckles on each shoulder strap and secure the buckles. Tighten the chest strap as necessary. **NOTE:** Ensure the orthosis is in a midline position with the spine and that it has not migrated upwards or downwards.

PANEL MODIFICATION FOR BRACE CUSTOMIZATION

NOTE: See customization worksheet for further details and documentation of brace customization for patient.

- The posterior panel of the brace can be modified using a regulated heat gun. Position the posterior panel in the appropriate lordotic curvature, apply heat, and allow to cool in that position. Modification is only recommended above the B and D ports as indicated on the inside of the posterior panel. The adhered hook may be removed during modification and reapplied once material has cooled. Let panel cool completely before applying brace to patient.

CARE INSTRUCTIONS

Remove panels from the brace and clean all soft goods by hand washing with mild soap and let air dry completely prior to reapplication. Anterior and posterior panels can be cleaned by wiping down with a disinfectant wipe.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DEROYAL® PRODUCTS ARE WARRANTED FOR NINETY (90) DAYS FROM THE DATE OF SHIPMENT FROM DEROYAL AS TO PRODUCT QUALITY AND WORKMANSHIP. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

PART#0-1991 | REVISED 3/2023

©2023 DEROYAL INDUSTRIES, INC.

All Rights Reserved. MAXAlign, DeRoyal, and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

Boa is a registered trademark of Boa Technologies Inc.

MAXALIGN® TLSO-D SPINE BRACE

	SINGLE PATIENT USE
	NON-STERILE
	NOT MADE WITH NATURAL RUBBER LATEX
RX ONLY	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings, and precautions before use. Correct application is essential for proper functioning of product.

INTENDED USE

The DeRoyal® MAXAlign® TLSO-D Spine Brace powered by the Boa® Fit System is intended to provide support and immobilization of the spinal column in the sagittal plane. The device is also intended to aid in relief and recovery following spinal surgery or an injury.

INDICATIONS

Post-operative thoracic/lumbar fusion, laminectomy, discectomy, compression fractures, degenerative disc disease, osteoporosis, single column instability, and facet syndrome.

CAUTIONS

- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.
- Inspect product for damage and/or missing components before use.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.**
- Prolonged usage can result in muscle atrophy in the user.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace.

INSTRUCTIONS FOR USE

SIZING

NOTE: Brace sizing is determined by circumferential waist measurement taken at the navel.

	ZONE 1	ZONE 2	ZONE 3	ZONE 4	ZONE 5	ZONE 6
MEASUREMENT	24" - 30"	30" - 36"	36" - 42"	42" - 48"	48" - 54"	54" - 60"

NOTE: For best results and ease of adjustment, brace sizing should be conducted on a large bed or table.

BRACE SIZING ADJUSTMENT:

- Pull out both Boa® Fit System dials to release lacing and release the front hook and loop closure to open the brace belt completely.
- Place the back of the brace flat with the outside of the

back panel facing upward. Ensure the left and right belt arms are spread outward evenly and are non-impinged.




- Unfasten the back-panel hook closures. Flip over the panel and belt arms to reveal the belt sizing adjustment mechanism.
- Pull up on the two (2) belt security hook flaps located on each side of the adjustment mechanism until they reach a 45-degree angle upward from the belt.
- Place one hand under the Belt Pull Strap near the hand indicator and use the other hand to pull up on the Belt Pull Strap until the appropriate size zone is reached on the left and right side of the belt. **NOTE:** Ensure that the sizing zones are even on both sides. Reapply the belt security hook closures to secure the belt.
- Align the Belt Pull Strap with the middle of the belt adjustment mechanism panel, ensure excess belt length is even and press down firmly. Fold the belt and back-flap over and reattach the back-flap hook closures at the top of the back panel.

DORSAL LUMBER EXTENSION (DLE) HEIGHT ADJUSTMENT:

- To adjust the Dorsal Lumbar Extension (DLE) height, unfasten the back-panel hook closures and flip over the back-panel and belt arms to access the DLE.
- Adjust the DLE higher or lower as necessary to fit patient. Position the DLE so that the yoke of the shoulder straps aligns with C7/T1 of patient.
- When appropriate height adjustment is complete, anchor the DLE using the hook sections placed in the middle and the top of the posterior panel.
- Fold the back-panel and belt over and reattach the back-flap hook closures at the top of the back panel.
- Detach the hook and loop anchors of the lower portion of the shoulder straps. These are positioned at the bottom of the posterior panel (Zone C).

BRACE APPLICATION:

- To apply the brace in the supine position, roll the patient on their side. Tuck the side closure of the belt behind the posterior panel and slide the closure under the patient's lower back so the panels are accessible and can be retrieved on the opposite side. Ensure the DLE straps are placed above the patient's shoulders so they can be

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