LIMB HOLDER

SUJECIÓN PARA EXTREMIDADES

ATTACHE-MEMBRES

GLIEDMASSENHALTER

REGGI-ARTO

SUPORTE DE MEMBROS

LEDEMAATHOUDER

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MANUFACTURED BY:
 DeRoyal Industries, Inc.
 200 DeBusk Lane
 Powell, TN 37849 USA
 888.938.7828 or (001) 865.938.7828
 www.deroyal.com

MedEnvoy Global B.V.
 Prinses Margrietplantsoen 33 – Suite 123
 2595 AM the Hague
 The Netherlands

UK Responsible Person
 MedEnvoy UK Limited
 85, Great Portland Street, First Floor
 London, W1W 7LT
 United Kingdom

PART#0-717G I REVISED 10/23

©2023 DeRoyal Industries, Inc.

All Rights Reserved. DeRoyal and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

EN: ENGLISH

| | |
|--|---|
| | SINGLE PATIENT USE |
| | NON-STERILE |
| | MEDICAL DEVICE |
| | NOT MADE WITH NATURAL RUBBER LATEX |
| | FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER. |

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

DeRoyal® Limb Holders are intended to be prescribed by a properly licensed practitioner to protect the patient from harming themselves, other patients, family members and/or staff.

CONTRAINDICATIONS

- Never use on a patient if an I.V. or wound site could be compromised by the device.
- DO NOT** use on a patient who is or may become highly aggressive or agitated.

⚠️ WARNINGS

- Keep straps and buckles out of patient's reach. Patient should not be able to remove device. Improper application or use of any restraint may result in serious injury or death.
- Inspect product for damage, missing components, or contamination before use. Ensure that all hook and loop, and buckle closures function properly.
- Use hospital approved bed safety guidelines. If securing the device to a bed always secure straps to the portion of the bed that moves with the patient. Do not secure straps to bed side rails.
- Check the patient regularly per hospital's policies and procedures to ensure circulation is not compromised.

CAUTIONS

- This product is to be fitted initially by a physician (or properly licensed practitioner) who is familiar with the purpose for which it is intended. The physician or practitioner is responsible for providing wearing instructions and cautions to other healthcare practitioners or healthcare providers involved in the patient's care.
- Consult physician (or properly licensed practitioner) immediately if patient experiences sensation changes, unusual reactions, swelling or increased pain while using this device. Be sure to follow your facility's policies and guidelines for frequency of patient monitoring.

INSTRUCTION FOR USE

SECURING THE LIMB HOLDER TO THE PATIENT

- Wrap the limb holder around the patient's wrist or ankle, ensuring that the cuff closure is on the outside of the device.
- Secure the hook and loop closure and if applicable, the buckle closure. **NOTE:** Two fingers should fit between the device and patient to ensure circulation is not compromised.
- Secure the strap(s) using quick-release ties or quick-release buckles depending on the style of the limb holder. Please refer to the step-by-step instructions for Quick-Release Ties.

SECURING THE STRAPS OF THE LIMB HOLDER TO THE BED:

- For limb holders with one strap, the strap should be attached to a single location on the portion of the bed that moves with the patient. Refer to quick-release ties and/or quick-release buckle instructions listed below.
- For limb holders with two straps, the two-point attachment option is available. One strap can be attached near the head of the bed and the other strap near the foot of the bed. The straps should be attached to the portion of the bed that moves with the patient. Refer to quick-release ties and/or quick-release buckle instructions listed below.
- Secure excess strap away from patient's reach.

QUICK-RELEASE TIES

Wrap the strap(s) around the anchor once. Take the loose end of the strap(s) and create a loop. Cross the loop over the other strap(s) that is wrapped around the anchor and feed it through opening as if tying a shoelace. Hold that loop and create another loop with the loose strap(s) and feed it through the first loop. Pull the final loop and the strap(s) attached to the cuff in opposite directions to form a quick-release tie.

QUICK-RELEASE BUCKLES

Release the female end of the buckle and wrap its strap around the anchor. Pull the female portion through the loop and connect the male and female ends of the buckle. Pull the excess strap to tighten.

CLEANING AND/OR MAINTENANCE

Please refer to hospital approved guidelines for cleaning and disinfecting of the device.

STORAGE AND TRANSPORT CONDITIONS

| | |
|--|--------------------------------|
| | KEEP DRY |
| | KEEP AWAY FROM SUNLIGHT |

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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