CATHETER STRAP

CORREA PARA CATÉTER

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KATHETERRIEMEN

FASCIA FISSACATETERE

CORREIA DE CATETER

KATHERBAND

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KATETER KAYIŞI

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CATHETER STRAP

	SINGLE PATIENT USE
	NON-STERILE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION

Please read all instructions, warnings, and cautions before use. Correct application is essential for proper functioning of the product. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INTENDED USE

The DeRoyal® Catheter Strap is intended to help secure catheter tubes in place.

WARNINGS

- Inspect product for damaged or missing components before use.
- Inspect the fit of the product on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
- Inspect the fit of the product around the catheter to ensure the catheter is adequately supported.
- Ensure circulation of fluid through the catheter is not compromised.
- Discontinue use and consult your prescriber if the product or its components break or become damaged, or cleaning is no longer sufficient to sanitize product.
- Do not use on a patient whose wound site could be compromised by the device.

CAUTIONS

- A physician or properly licensed practitioner (a “prescriber”) who is familiar with the use and purpose of this product must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this product to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.
- Consult prescriber immediately if patient experiences sensation changes, unusual reactions, swelling or increased pain while using this product. Discontinue use if pressure injuries develop.
- Take special care if the user is diabetic or has poor circulation. These users may have decreased skin sensitivity and are at greater risk for poor peripheral circulation and pressure injuries.

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INSTRUCTION FOR USE

1. Measure the circumference of the patient's thigh to determine the appropriate size.
2. Position the catheter strap behind the patient's thigh with the wide blue tab on the medial side of the patient's leg.
3. Wrap the elastic catheter strap around the thigh and secure with closure.
4. Position the catheter vertically atop the sewn end of the blue tab.
5. Wrap the tab around the catheter strap and secure the closure.
6. Place the small white foam piece on top of the exposed tab that wraps around the catheter strap.
7. The catheter should have the ability to move with the patient to help prevent bladder traction.

SIZE	DIMENSIONS
Universal	22" L x 2" W (56 x 5 cm)
Bariatric	38" L x 2" W (100 x 5 cm)

CLEANING AND/OR MAINTENANCE

Refer to hospital protocols and guidelines.

STORAGE AND TRANSPORT CONDITIONS

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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All Rights Reserved. DeRoyal, the DeRoyal logo, and Improving Care. Improving Business. are trademarks or registered trademarks of DeRoyal Industries, Inc.