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STATIC-PRO[®] WRIST
STATIC PROGRESSIVE ORTHOSIS

SINGLE PATIENT USE
NON-Sterile
NOT MADE WITH NATURAL RUBBER LATEX
FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.
RX ONLY

IMPORTANT INFORMATION

Please read all instructions, warnings and cautions before use. Correct application is essential for proper product function and to reduce the risk of injury or re-injury inherent with the use of any brace. Use only on the person it was provided to by a healthcare professional and only for the use it was intended.

INDICATIONS

Not to be used if any of the following conditions exist: Severe spasticity, severe osteoporosis, and/or thromboembolitis.

CAUTIONS

- A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.
- Use only as directed by medical professional.
- Prolonged usage can result in muscle atrophy in the user. The user may require (1) a weaning period from the brace to build muscle strength (2) regular physical therapy exercises to maintain the strength of the wrist (3) a limitation on duration of use or (4) some combination of these.
- Inspect brace for damaged or missing components before use.
- Discontinue use and consult your prescriber if the brace or its components break or become damaged.
- Inspect the fit of the brace on the patient to ensure fit is proper. **DO NOT OVERTIGHTEN.** Check regularly to ensure circulation is not compromised.
- Consult your prescriber immediately if you experience sensation changes, unusual reactions, swelling or increased pain while using this brace. Discontinue use if pressure injuries develop.
- Take special care if the user is diabetic or has poor circulation as these users may have decreased skin sensitivity and are at greater risk for poor peripheral circulation and pressure injuries.

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
INSTRUCTIONS FOR USE

IMPORTANT: Before applying the orthosis, disengage the gearing

- ENGAGE GEARING MECHANISM:** Pull the black engagement knob (1) outward and pivot the position adjustment dial (2) forward towards the hand. The black engagement knob will click into place. **NOTE:** To fully engage the gearing mechanism, the position adjustment dial (2) may need to be rotated slightly.


- ULNAR/RADIAL DEVIATION:** Adjust for ulnar/radial deviation by using the hex key wrench, and loosening the two screws (10) on the ulnar/radial plate (9). Adjust to desired position, and re-tighten the screws.

- ADJUST THE ORTHOSIS POSITION:** Rotate the position adjustment dial (2) until the joint's soft tissues are stretched to patient tolerance. **NOTE:** See Patient Protocol Card to determine the appropriate amount of stress, and wearing time.

- REMOVING ORTHOSIS:** Disengage the gearing mechanism by pulling out the engagement knob (1), and at the same time pivot the position adjustment dial (2) toward the elbow. The handcuff should move freely. **NOTE:** The adjustment dial should be positioned below the wrist in order for the device to achieve full extension.

- SIZING THE FOREARM CUFF:** The forearm cuff (3) is designed to be trimmed to fit each patient. Unfasten forearm straps (4), and using scissors, trim from both sides of the cuff so that the ends of the cuff do not overlap when on the forearm.

- DONNING THE UNIT:** Slide the hand and forearm into the orthosis, with the mechanism on the medial (inside) of the forearm. Position the axis of the hinge centered at the wrist joint. Mold the forearm cuff (3) around the forearm, and fasten the two straps with the hook and loop closure for a snug, comfortable fit. **NOTE:** Pulling the welded end of each strap will

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY

DeRoyal[®] products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROVAL'S WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

STORAGE AND TRANSPORT CONDITIONS

KEEP DRY
KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

CAUTIONS

• A physician or properly licensed practitioner (a "prescriber") who is familiar with the use and purpose of this brace must fit it to the user. The prescriber has a duty to provide wearing instructions and risks related to the use of this brace to other healthcare practitioners treating the users and the users themselves, including duration of use. The instructions provided in this sheet do not supersede hospital protocol or direct orders of the prescriber.

• Use only as directed by medical professional.

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All Rights Reserved. DeRoyal, the DeRoyal logo,
and Static-Pro are trademarks or registered
trademarks of DeRoyal Industries, Inc.