
DeRoyal uses the symbols and meanings from standards ISO 15523-1. These symbols are placed next to the text explaining their meaning in this instructions for use (IFU), a complete symbols glossary is available online at deroyal.com/symbols or by contacting customer service.

EN: FOLEY CATHETER WITH TEMPERATURE SENSOR

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PT: CATETER DE FOLEY COM SENSOR DE TEMPERATURA

Manufactured by:
DeRoyal Industries, Inc.
200 DeBusk Lane
Powell, TN 37849 USA
888.938.7828 or (001) 865.938.7828
www.deroyal.com

PART#74-14903
Revised 2/2025
©2025 DeRoyal Industries, Inc.
All Rights Reserved. DeRoyal and the DeRoyal logo are trademarks or registered trademarks of DeRoyal Industries, Inc.

EN: ENGLISH

FOLEY CATHETER WITH TEMPERATURE SENSOR

	MR CONDITIONAL
	DO NOT RESTERILIZE
	DO NOT REUSE
	STERILIZED WITH ETHYLENE OXIDE
	NOT MADE WITH NATURAL RUBBER LATEX
	FEDERAL U.S.A. LAW RESTRICTS THIS DEVICE TO SALE OR USE BY OR ON THE ORDER OF A PHYSICIAN OR PROPERLY LICENSED PRACTITIONER.

IMPORTANT INFORMATION
Please read all warnings and instructions before use. Correct application is essential for proper product function and to reduce the risk of injury.

INTENDED USE
The DeRoyal® Foley Catheter with Temperature Sensor is to be used for drainage of the urinary bladder and simultaneous monitoring of core body temperature.

INTENDED USERS
The DeRoyal® Foley Catheter with Temperature Sensor is intended to be used only by a licensed practitioner or trained healthcare professionals.

PATIENT TARGET GROUP(S)
The intended patient population is adults and pediatric.

DESCRIPTION
The Foley Catheter with Temperature Sensor consists of a temperature-sensing wire set secured inside a silicone catheter. The catheter's balloon inflation system is terminated with a luer tip syringe-actuated, self-sealing valve. The device is available in 8- 18 French size. The French size and balloon volume are designated on the catheter.

The wire set is secured within the lumen of the foley catheter. The temperature sensor is located near the distal tip of the catheter. The connector at the proximal end interfaces with DeRoyal-approved cables for connection to a monitoring unit. The temperature monitoring sensor is available in 400 series configurations for connection to a compatible monitoring unit. Consult the user manual for the patient monitor to ensure compatibility.

INDICATIONS
Foley catheterization is recommended when continuous bladder drainage is required in the absence of voluntary micturition, after surgery or trauma of the pelvic organs, in the presence of obstruction or paralysis, and as a means to evaluate urine output in patients having fluid management problems.

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CONTRAINDICACIONES
None known.

! WARNING

- The patient should be routinely monitored for urine output to minimize risk of urine retention.
- Use aseptic technique in patient preparation and catheter care.
- Follow accepted medical techniques during removal of catheter to minimize risk for patient infection or injury.
- Carefully inspect the packaging and device for damage or defect. **DO NOT** use if the packaging is open or the device is damaged.
- DO NOT** use device on a patient with weak bladder walls.
- This device has not been evaluated for reprocessing or re-sterilization. Reprocessing and/or re-sterilization may damage the device, rendering it unusable and/or may lead to device failure.
- DO NOT** cut, damage, or augment this device in any way. Doing so may lead to injury to the patient/clinician and/or may damage the internal structures of the probe causing leaks and loss of function.
- Chronic catheter use increases the risk of bladder retention, urinary tract infection or other adverse event. Minimize catheterization duration following accepted clinical practice. If resistance to flow in the catheter develops, then the catheter should be changed.
- DO NOT** store device in extreme temperature or a moist environment.
- This device is for single patient use only. It is not intended for disinfection and subsequent re-use, which may result in device failure and/or create the risk of contamination.
- DO NOT** use an expired device.
- Prior to removing the catheter, deflate the balloon by inserting the luer tip of an empty syringe into the valve and aspirating the solution completely. If the balloon does not collapse after attempting this step, then sever the valve arm. If balloon still does not deflate, call physician.

ADVERSE REACTIONS
While adverse reactions are rare, the following are potential adverse reactions associated with the use of temperature probes during insertion or while the device is in use:

- allergic reaction to materials
- bladder spasms
- electric burns due to aberrant electro-surgical current
- encrustations with urinary salts
- kidney or bladder stones
- septicemia
- urethral discharge
- urethritis
- urethrocutaneous fistulae
- urinary tract infection or irritation

CAUTIONS

- The clinician is responsible for determining the proper product size and the patient's suitability for the safe and effective use of this product.
- Take proper care to not damage the insulation prior to or during use.
- If temperature does not read properly, check all connections. If

problem persists, discard and use a new probe. Probe may take some time to stabilize. If the response time is excessive, discard and replace.

- Ensure the cable and connector do not get wet.
- Prior to insertion, test the balloon inflation and deflation. The balloon should be fully deflated before insertion.
- Follow the recommended balloon inflation volume printed on the product and use sterile water for inflation.
- Select the appropriate size catheter using expert clinical and medical judgment.
- Lubricate catheter with a suitable water soluble lubricant prior to insertion using accepted medical techniques.
- DO NOT** use a needle to puncture the catheter wall for aspiration. Leaks may occur.
- DO NOT** intertwine the cables, especially the monitor cables with the electro-surgical unit's cables.
- The operation of the patient temperature monitor may be temporarily affected during electro-surgical activations. Unusual temperature readings should be checked.
- The catheter is intended to be used with DeRoyal-approved interface cables. Use with an incompatible cable may affect performance.

DIRECTIONS FOR USE

- Use aseptic technique in patient preparation and catheter care.
- Verify the compatibility of the catheter, interface cable, and patient monitor before use.
- Remove the Foley Catheter with Temperature Sensor from its sterile package using aseptic technique.
- Test the balloon inflation and deflation by inserting a luer tip syringe filled with the recommended amount of sterile water into the valve. Be sure to fully deflate the balloon before insertion.
- Lubricate the catheter with a suitable water soluble lubricant and insert using acceptable medical techniques.
- Make sure that the catheter is draining properly before inflating the balloon.
- To inflate the balloon, insert a luer tip syringe filled with the recommended amount of sterile water into the valve. Depress the plunger. The valve will seal when the syringe tip is withdrawn.
- Attach drainage tubing and bag.
- Secure the catheter, drainage tubing and bag in place by following accepted medical technique.
- Connect the temperature probe to the appropriate reusable cable from DeRoyal. Connect the cable to the patient monitor. Secure the cable using the drape clip.
- Follow the directions for use of the patient temperature monitor instrument.
- Prior to removing the catheter, deflate the balloon by inserting the luer tip of an empty syringe into the valve and aspirating the solution completely.
- If balloon does not collapse after attempting to perform #12 above, sever the valve arm. If balloon still does not deflate, call physician.
- Remove the catheter following accepted medical techniques.
- Disconnect the probe (at the connector) and discard probe. To disconnect, grasp both connectors firmly and pull. **DO NOT** pull on cable or wire.

MR COMPATIBILITY
The Foley Catheter with Temperature Sensor was determined to be MR Conditional according to the terminology specified in the American Society for Testing and Materials ("ASTM") International

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designation: Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment, and other related standards.

This Foley Catheter with Temperature Sensor is safe for conditional use in MRI imaging when utilized in compliance with the instructions located here: www.deroyal.com/foley-catheter-mri-compatibility or by calling: 1-800-251-9864. **WARNING: USE OUTSIDE OF THESE PARAMETERS COULD RESULT IN SERIOUS INJURY OR DEATH.**

OPERATING SPECIFICATIONS

- Operating Temperature: 25° C to 45° C
- Rated Output Range: 35° C to 42° C
- Rated Extended Output Range: 25° C to 35° C and 42° C to 45° C
- The reference body site is the core body temperature.
- The measuring site is the Bladder.
- Accuracy in Rated Output Range: +/- 0.2° C
- Accuracy in Rated Extended Output Range: +/- 0.3° C
- Time response:
 - The temperature probe has a direct mode of operation.
 - The maximum heating transient time is 1 minute
 - The maximum cooling transient time is 30 seconds

DISPOSAL
Dispose of any used device according to local, state, and federal laws and regulations. For safe disposal of devices, follow your facility's protocol.

	KEEP DRY
	KEEP AWAY FROM SUNLIGHT

In addition to the competent authority in the country where the patient resides, serious incidents must be reported to DeRoyal Industries, Inc.

WARRANTY
DeRoyal® products are warranted for one hundred twenty (120) days from the date of shipment from DeRoyal as to product quality and workmanship. **DEROYALS WRITTEN WARRANTIES ARE GIVEN IN LIEU OF ANY IMPLIED WARRANTIES, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.**

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