




## Non-Sterile Containers & Trays

for OEM Manufacturers & Kit Packers



*A message from all of us at Hospital Disposables...*

*Our goal is to become your partner of choice for disposable plastic medical containers and trays. We are committed to maintaining operational excellence that will add value to your company.*

*\* All products packed non-sterile.*

## CONTENT

Product Listing	Page
Sponge Bowls . . . . .	3
Guidewire Bowl . . . . .	3
Ring Basins . . . . .	4
Basin Strainers . . . . .	4
Emesis Basin . . . . .	5
Placenta Basin . . . . .	5
Rectangular Basin . . . . .	5
Pitcher . . . . .	6
Medicine Cups . . . . .	6
Trays . . . . .	7

## HOSPITAL DISPOSABLES

### Sponge Bowls

- Designed to hold saline solutions, sponges, and other medical supplies
- Graduated in oz, cc, and ml for easy measurements
- Optional lid
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-0250	Bowl, Denture Cup, 8 oz (250 cc)	600/Cs
395-9025	Lid, Denture Cup, 8 oz (250 cc)	300/Cs
395-0500	Bowl, 16 oz (500 cc)	390/Cs
395-9050	Lid, 16 oz (500 cc)	270/Cs
395-1000	Bowl, 32 oz (1000 cc)	160/Cs
395-9100	Lid, 32 oz (1000 cc)	160/Cs
395-1500	Bowl, 48 oz (1500 cc)	160/Cs
395-9150	Lid, 48 oz (1500 cc)	160/Cs
395-2500	Bowl, 80 oz (2500 cc)	160/Cs
395-9250	Lid, 80 oz (2500 cc)	120/Cs




### Guidewire Bowl

- Built-in lip and recessed bottom to hold guidewire in place and allow for easy removal
- Graduated in oz, cc, and ml for easy measurements
- Optional Lid
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-GW2000	Bowl, 80 oz (2500 cc)	32/Cs
395-GW2000CT	Bowl, 80 oz (2500 cc)	160/Cs
395-9250	Lid, 80 oz (2500 cc)	120/Cs

\* Patent No. D433,130



# Non-Sterile Containers & Trays

## HOSPITAL DISPOSABLES

### Ring Basins

- Designed to hold saline solutions, instruments, and other medical supplies
- Fits all standard ring stands
- Graduated in qt, cc, and ml for easy measurements
- ETO and gamma sterilizable unless otherwise noted
- Not made with natural rubber latex

Product No.	Description	Quantity
395-6000	Basin, 6 qt (6000 cc)	60/Cs
395-7000*	Basin, 7 qt, (6500 cc)	38/Cs

\* **ETO Sterilizable only – contact Hospital Disposables for information on gamma sterilizable product.**  
**Graduated in cc only.**



### Basin Strainers

- Facilitates removal of contaminated instruments from basins
- Does not interfere with normal ring basin use
- Fits standard ring basins
- Graduated in qt, cc, and ml for easy measurements
- Optional “shower cap” lid
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-ST6000	Strainer, Blue, 6 qt (6000 cc)	60/Cs
395-ST6000-RD	Strainer, Red, 6 qt (6000 cc)	60/Cs
395-SB9600	Lid, 6 qt (6000 cc)	100/Cs



## HOSPITAL DISPOSABLES

### Emesis Basin

- Flexible design with rounded edges
- Graduated in oz, cc, and ml for easy measurements
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-0700	Basin, 700 cc	320/Cs



### Placenta Basin

- Graduated in oz, cc, and ml for easy measurements
- Optional lid
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-3000	Basin, 3000 cc	100/Cs
395-9300	Lid, 3000 cc	100/Cs



### Rectangular Basin

- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
395-8000	Basin, 8 qt (8000 cc)	50/Cs



# Non-Sterile Containers & Trays

## HOSPITAL DISPOSABLES

### Pitcher

- Clear material for easy viewing of contents
- Pour spout for easy transfer of liquids
- Graduated on both sides in oz and cc for easy use by right or left-handed person
- Handle for easy use
- Not made with natural rubber latex

Product No.	Description	Quantity
395-1200*	Pitcher, 40 oz (1200 cc)	70/Cs
395-1200CT*	Pitcher, 40 oz (1200 cc)	1344/Cs

\* ETO Sterilizable only



### Medicine Cups

- Tapered design
- Graduated in oz and cc
- Marked on two sides every ¼ oz and 5 cc
- Optional lid
- Not made with natural rubber latex

Product No.	Description	Quantity
395-0060*	Cup, Clear, 2 oz (60 cc)	1380/Cs
395-0060-BL*	Cup, Blue, 2 oz (60 cc)	1380/Cs
395-0060-GN*	Cup, Green, 2 oz (60 cc)	1380/Cs
395-0060RD*	Cup, Red, 2 oz (60 cc)	1380/Cs
395-0060-YL*	Cup, Yellow, 2 oz (60 cc)	1380/Cs
395-9006*	Lid, 2 oz (60 cc)	1380/Cs

\* ETO Sterilizable only




## HOSPITAL DISPOSABLES

### Trays

- Injection molded
- Smooth edges that maintain the integrity of wraps, pouches, and bags
- Thick corners that do not crack or break
- “Clean” product that does not need wiping
- ETO and gamma sterilizable
- Not made with natural rubber latex

Product No.	Description	Quantity
396-0542	Small Deep, 4¾" x 5½" x 2"	600/Cs
396-0952	Large Deep, 5½" x 9¾" x 2"	300/Cs
396-0541	Small Shallow, 4¾" x 5½" x 1"	600/Cs
396-0951	Large Shallow, 5½" x 9¾" x 1"	300/Cs
396-0904	2 Compartment Prep, 4½" x 9½" x 1¾"	300/Cs
396-0910	3 Compartment Prep, 9½" x 10" x 2"	150/Cs
396-0911	J-Wire, 9¾" x 11" x 2"	150/Cs
396-1219	Small Mayo, 12¾" x 19¼"	40/Cs
396-1622	Large Mayo, 16½" x 21¼"	30/Cs
396-1115	Small Platform, 10¾" x 14¾" x 1"	80/Cs
396-1416	Large Platform, 14½" x 16¼" x 1"	50/Cs





## Hospital Disposables

A Division of DeRoyal Industries, Inc.

200 DeBusk Lane  
Powell, TN 37849  
800-327-0450  
865-938-7828  
E-mail: [oeminfo@deroyal.com](mailto:oeminfo@deroyal.com)

©2019 DeRoyal | Reprint # 0-1534 | Rev. 04/19