

TEMPERATURE THERAPY

UNITS & BLANKETS





JetStream® Hot/Cold Temperature Therapy

- Patented ATC™ measures and regulates temperature to provide therapeutic hot or cold therapy without manual adjustment
- Intuitive LED light lets user know when in therapeutic zone and when to add ice or hot water
- Portable system suitable for any medical setting or treatment at home
- Lightweight, low-profile design and Quick Start Guide are patient friendly
- Both hot and cold therapy in one unit
- Therapy Blankets optimize temperature distribution, insulation and comfort
- Blanket Wraps increase insulation and ensure precise placement of blankets

T700 Combo Units	Non-Sterile	Sterile
Hip/Lumbar	T725NS	T725ST
Knee/Shoulder	T735NS	T735ST
Knee/Shoulder (w/attached straps)	T736NS	T736ST
Shoulder/Knee (w/attached straps)	T737NS	T737ST

COOLJET™



CoolJet™ **Cold** Temperature Therapy

- Patented ATC™ measures and regulates temperature to provide therapeutic cold therapy without manual adjustment
- Intuitive LED light lets user know when in therapeutic zone and when to add ice
- Portable system suitable for any medical setting or treatment at home
- Therapy Blankets optimize temperature distribution, insulation and comfort
- Blanket Wraps increase insulation and ensure precise placement of blankets
- Intuitive screw top opening

T800 Combo Units	Non-Sterile	Sterile
Hip/Lumbar	T825NS	T825ST
Knee/Shoulder	T835NS	T835ST
Knee/Shoulder (w/attached straps)	T836NS	T836ST
Shoulder/Knee (w/attached straps)	T837NS	T837ST



Blankets & Wraps

Blankets	Blanket Dimensions	Strap Lengths	Non-Sterile	Sterile
Arm/Leg	10" x 5"	-	T5015NS	T5015ST
Hip/Lumbar	10.5" x 11.5"	-	T5025NS	T5025ST
Knee/Shoulder	10.5" x 11"	24" detached (2) 36" detached (1)	T5035NS T5035NS-10 (10 per/case)	T5035ST T5035ST-10 (10 per/case)
Knee/Shoulder-Large	12.5" x 11"	24" detached (2) 36" detached (1)	T5035NS-L T5035NS-L-10 (10 per/case)	T5035ST-L T5035ST-L-10 (10 per/case)
Knee/Shoulder (w/ attached straps)	10.5" x 11"	24" attached (2)	T5036NS T5036NS-10 (10 per/case)	T5036ST T5036ST-10 (10 per/case)
Knee/Shoulder-Large (w/attached straps)	12.5" x 11"	24" attached (2)	T5036NS-L T5036NS-L-10 (10 per/case)	T5036ST-L T5036ST-L-10 (10 per/case)
Shoulder/Knee (w/ attached straps)	10.5" x 11"	36" attached (1) 12" attached (1)	T5037NS T5037NS-10 (10 per/case)	T5037ST T5037ST-10 (10 per/case)
Shoulder/Knee - Large (w/attached straps)	12.5" x 11"	36" attached (1) 12" attached (1)	T5037NS-L T5037NS-L-10 (10 per/case)	T5037ST-L T5037ST-L-10 (10 per/case)
Hand/Wrist	3.5" x 5"	-	T5055NS	T5055ST
Cervical	10" x 21"	-	T5075NS	T5075ST
Foot/Ankle	10" x 11.5"	-	T5085NS	T5085ST

Foam Wraps	Product #
Chest/Lumbar	T6100
Knee	T6110
Shoulder	T6120
Hips	T6160

DeRoyal®, DeRoyal logo, Improving Care. Improving Business.®, and JetStream® are registered trademarks of DeRoyal Industries, Inc.

CoolJet™ & ATC™ are trademarks of DeRoyal Industries, Inc.

Copyright© 2016, DeRoyal. All rights reserved.



Improving Care. Improving Business.™

