

*Solace*TM

Seek Comfort.

Seek Support.

Seek Solace.



LOW PROFILE BACK BRACE

 DeRoyal[®]

Solace™

LOW PROFILE BACK BRACE

DeRoyal's Solace Low Profile Back Brace combines quality, comfort, and a low profile design to provide your patients with the relief they need.

Features & Benefits

- Lace up closure system easily adjusts providing maximum compression
- Soft, breathable material provides a cool, comfortable fit increasing patient compliance
- Rigid anterior and posterior plastic inserts provide support to the intra-abdominal and viscera region
- Low profile, lightweight design fits easily under patient's clothing
- Four sizes available to ensure a comfortable fit
- Plastic D-rings for ease of application and removal

Indications

- Chronic low back pain
- Lumbar muscle weakness
- Lumbar strain & sprain
- Mechanical or discogenic lumbar pain

Sizing Information

Product #	Size	Waist Circumference
13400005	Small	26" - 34"
13400006	Medium	32" - 40"
13400007	Large	38" - 46"
13400008	X-Large	44" - 51"

Front



Side



Back



Diverse. Dynamic.

