

Cellulose/Polyester Non-Woven Wiper



DERTEX™

Will not scratch surfaces. Use either in wet or dry applications. Class 100.

DERTEX™ Cellulose/Polyester Wipers are General Purpose Non-Woven Wipers, made from a hydroentangled blend of 55% Cellulose and 45% Polyester Fibers. This fiber combination provides low particulate and pure wipe performance.

Material is soft, durable and highly absorbent.
Will not scratch surfaces.
Use either in wet or dry applications.

Economical and efficient general purpose wiper that offers strength, economy and cleanliness for a synthetic wiper.



6" x 6" Sheets NCP-668QFB

9" x 9" Sheets NCP-968

Poly/Cellulose Non-Woven Wiper

Item #	Description	Packaging
NCP-668QFB	6" x 6" (15 cm x 15 cm) Non-Woven, Cellulose/Polyester Blend Wipers 68 gsm Weight	100 Wipers/Box 10 Inner Bags/Case
NCP-968	9" x 9" (23 cm x 23 cm) Non-Woven, Cellulose/Polyester Blend Wipers 68 gsm Weight	300 Wipers/Bag 10 Inner Bags/Case
NCP-1268	12" x 12" (31 cm x 31 cm) Non-Woven, Cellulose/Polyester Blend Wipers 68 gsm Weight	150 Wipers/Bag 10 Inner Bags/Case
NCPC-660QF	6" x 6" (15 cm x 15 cm) CREPED, Qtr. Folded, Non-Woven Cellulose/Polyester 60 gsm Weight	50 Wipers/Bag 12 Inner Bags/Case

Applications / Features

- Designed for General Purpose wiping where economy and cleanliness are critical.
- Ideal for applications where high absorbency and low extractables are a major requirement.
- Suitable for cleaning critical components and lab equipment.



12" x 12" Sheets NCP-1268



Qtr. Fold 6" x 6" Sheets NCPC-660QF

DERTEX™

To Order Call: 800.845.1062, Fax: 803.424.5112 or E-mail: Textiles@DeRoyal.com ©2013 DeRoyal Inc.

