


Heel Suspension Boot



Designed to completely off-load pressure from the heel, DeRoyal's Heel Suspension Boot helps prevent heel pressure ulcers.

Features & Benefits

- Smooth, breathable foam construction
- Sloped wedge inside boot completely suspends the heel
- Ventilation holes help promote air circulation and dissipate heat
- Soft straps with hook closures for application and easy adjustability
- Posterior panel helps reduce friction against bed surface
- Extra pad attaches to side of boot to help prevent hip rotation
- Universal in size
- Latex safe

Indications

- Heel pressure ulcers
- Hip rotation
- Hip fractures (pre and post-operative)

Ordering Information

M60-300

*Packaged 6/Case



200 DeBusk Lane
Powell, TN 37949
www.deroyal.com
800.251.9864 • 865.938.7828
Reprint# LIT0-2097
rev. 02/10