



FULL SHELL 2

WALKER BOOT

COMFORT for the way you **LIVE**

CONTROL for the way you **PLAY**

CONFIDENCE in how you **RECOVER**

FULL SHELL Walker 2

FEATURES & BENEFITS

- Built in medial and lateral air pumps allow easy operation and instant air bladder control.
- Air bladders accommodate swelling and bulky dressings.
- Wide foot bed accommodates varying patient foot anatomies.
- Textured tread, non-skid sole provides maximum shock absorption and grip.
- Rocker bottom mimics natural gait.
- Unique adjustable toe guard for increased protection.
- Additional pads included for extra comfort where needed.



INDICATIONS: Moderate to severe ankle sprains and strains | Stable foot and/or ankle fractures | Syndesmosis sprain | Post ORIF | Achilles tendonitis

SIZE	MALE SHOE SIZE	FEMALE SHOE SIZE	STANDARD	SHORT
X-SMALL	up to 4	up to 5	15500041	15500046
SMALL	4 - 7	5 - 8	15500042	15500047
MEDIUM	7 - 10	8 - 11	15500043	15500048
LARGE	10 - 13	11 - 15	15500044	15500049
X-LARGE	13 +	15 +	15500045	15500050



Improving Care. Improving Business.®

DeRoyal Industries, Inc.
 200 DeBusk Ln | Powell, TN 37849
 888.938.7828 | www.deroyal.com
 Revision 6/14 | Reprint # 0-2347

