


Kittner

Laparoscopic Dissection Sponge

- Recommended for laparoscopic, endoscopic thoracoscopic (VATs) and robotic procedures
- Physician designed
- X-ray detectable
- Tapered edges
- High quality rayon fabric
- Patent pending
- Sterile
- Not made with natural rubber latex

The DeRoyal® Kittner Laparoscopic Dissection Sponge is an x-ray detectable sponge intended to be inserted through a trocar by a licensed practitioner for tissue dissection, mobilization, fluid absorption, moisture retention and tissue protection during various surgical procedures.



Ordering Information

30-305K | Kittner Laparoscopic Dissection Sponge | 5/pk 50/Cs

DeRoyal, the DeRoyal logo, Improving Care. Improving Business, are registered trademarks of DeRoyal Industries, Inc.
© 2017, DeRoyal. All rights reserved.

DeRoyal Industries, Inc.
200 DeBusk Lane | Powell, TN 37849
1.800.DEROYAL or 865.938.7828

www.deroyal.com | Rev. 12/18 Reprint #0-2444



Improving Care. Improving Business.®