



Suction Tubing & Handles

www.deroyal.com



Improving Care. Improving Business.®

Suction Tubing with Connector

- Molded female connectors provides secure connection
- Male connector included for coupling additional tubes
- Minimal coil memory allows for easy handling and drapability
- Color coded connectors for size identification
- Available in economical non-sterile 100ft. rolls
- Not made with natural rubber latex



3/16" Tubing with Green Female Connectors

Product #	Description	Quantity
71-C501*	3/16" x 1.5' (5 mm x .46 m)	100/Cs
71-C506*	3/16" x 6' (5 mm x 1.8 m)	50/Cs
71-C510*	3/16" x 10' (5 mm x 3 m)	30/Cs
71-C512*	3/16" x 12' (5 mm x 3.7 m)	30/Cs
71-C520*	3/16" x 20' (5 mm x 6.1 m)	25/Cs
71-C5100NS*	3/16" x 100' (5 mm x 30 m)	1 roll/Bx




1/4" Tubing with Light Blue Female Connectors

Product #	Description	Quantity
71-C606*	1/4" x 6' (6 mm x 1.8 m)	45/Cs
71-C610*	1/4" x 10' (6 mm x 3.0 m)	25/Cs
71-C612*	1/4" x 12' (6 mm x 3.7 m)	25/Cs
71-C620*	1/4" x 20' (6 mm x 6.1 m)	20/Cs
71-C6100NS*	1/4" x 100' (6 mm x 30 m)	1 roll/Bx



9/32" Tubing with Dark Blue Female Connectors

Product #	Description	Quantity
71-C706*	9/32" x 6' (7 mm x 1.8 m)	40/Cs
71-C710*	9/32" x 10' (7 mm x 3.0 m)	20/Cs
71-C712*	9/32" x 12' (7 mm x 3.7 m)	20/Cs
71-C720*	9/32" x 20' (7 mm x 6.1 m)	15/Cs
71-C7100NS*	9/32" x 100' (7 mm x 30 m)	1 roll/Bx



⚠️ WARNING - This product can expose you to chemicals including di(2-ethylhexyl)phthalate, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

Suction Handles



Flexible Suction Handles

- Constructed of a malleable non-breakable vinyl with smooth inner lumen
- Highly flexible for precise suctioning
- Transparent for visual monitoring and inspection during use
- Available in three tip styles for different volumes of suction (fine, regular, high)
- Side eyes help to reduce the chance of tissue grab and/or trauma
- Vented option available with regular capacity
- Not made from natural rubber latex



Flexible Suction Handle Tip

Product #	Description	Quantity
71-D860	Fine capacity	50/Cs
71-D861	Regular capacity,	50/Cs
71-D861V	Regular capacity, vented	50/Cs
71-D862	High capacity	50/Cs

Poole Suction Handles



Poole with Outer Sheath




Poole Handle with Tubing

- Constructed from a shatter-resistant plastic
- Designed to remove pooled blood and fluids during surgical procedures
- Features three multipurpose positions that convert from sump action to pinpoint capabilities
- Outer sheath has 132 large eyes for rapid suctioning while reducing tissue grab
- Inner cannula provides pinpoint suctioning
- Universal connector for easy connection to a variety of tubing sizes
- Available with pre-attached tubing
- Not made from natural rubber latex

Product #	Description	Quantity
71-D8770	Poole Suction Handle	50/Cs
71-D8776*	Poole Suction Handle with 1/4" (6mm) x 6' (1.8m) tubing	20/Cs

Frazier Suction Handles



Frazier Suction Handle Tip

- Ideal for surgical procedures that require fine and delicate suctioning
- Malleable shaft can be bent to desired position without occlusion
- Obturator can be used to facilitate shaping of shaft and/or removal of obstructions that may become trapped during suctioning
- Obturator must be removed before suctioning begins
- Control vent allows for intermittent or continuous suction control
- Universal connector for easy connection to a variety of tubing sizes
- Not made from natural rubber latex


Product #	Description	Quantity
71-D871	8 French	50/Cs
71-D872	10 French	50/Cs
71-D873	12 French	50/Cs

WARNING - This product can expose you to chemicals including di(2-ethylhexyl)phthalate, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

Yankauer Suction Handles

- Constructed of a shatter-resistant plastic
- Smooth inner surface for maximum flow
- Slip-resistant and lightweight, provides optimal control and maneuverability
- Translucent material allows for visual monitoring and inspection during use
- Available in three tip styles: bulbous, tapered bulbous and straight
- Bulbous tip is rounded to help minimize tissue grab and/or trauma
- Tapered bulbous tip helps minimize tissue grab and/or trauma, while allowing for better visibility
- Straight tip provides less bulk, allowing for better visibility
- Each tip has side eyes to help reduce tissue grab
- Available with control vent for intermittent or continuous suction control
- Available with pre-attached tubing
- Not made from natural rubber latex

Product #	Description	Quantity
71-D886	Bulbous tip, non-vented	50/Cs
71-D886V	Bulbous tip, vented	50/Cs
71-D880	Tapered bulbous tip, non-vented	50/Cs
71-D882	Tapered bulbous tip, vented	50/Cs
71-D887	Straight tip, non-vented	50/Cs
71-D887V	Straight tip, vented	50/Cs
71-D883A*	Tapered bulbous tip, non-vented, with ¼" (6mm) x 6' (1.8m) tubing	50/Cs
71-D884*	Tapered bulbous tip, non-vented, with ¼" (6mm) x 9' (2.7m) tubing	25/Cs
71-D885*	Tapered bulbous tip, non-vented, with ¼" (6mm) x 12' (3.7m) tubing	25/Cs



WARNING - This product can expose you to chemicals including di(2-ethylhexyl)phthalate, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

DeRoyal, the DeRoyal logo, and Improving Care. Improving Business.
are registered trademarks of DeRoyal Industries, Inc.
PRUventor and Silbac are trademarks of DeRoyal Industries, Inc.
©2019, DeRoyal. All rights reserved.

DeRoyal Industries, Inc. | 200 DeBusk Lane, Powell, TN 37849 USA
888.938.7828 or 865.938.7828 | Customer Service: 800.251.9864 | www.deroyal.com
Reprint # 0-2448 | Rev. 01/19

