


Enhanced Foot-Drop Support



Features and Benefits

- **Static Bilateral Straps** – Keeps the foot in neutral position
- **Structured Foot Rest** – Additional foot drop support
- **Foot Securement Strap (optional)** - Improves patient compliance
- **Silbac™ Coating** - Anti-Microbial technology
- **Leg Securement Strap** – Adjustable and customizable
- **Anti-Rotation Wedge** – Helps prevent foot rotation
- **SCDI/ICD Exits:** Outlined by white stitching

Clinically Proven to Help Prevent Heel Pressure Injuries*



*data on file

⚠️ *WARNING: This product can expose you to chemicals including N,N-Dimethylformamide, which is known to the State of California to cause cancer. For more information, go to www.P65Warnings.ca.gov.

Ordering Information

Product #	Description	Length	Quantity
M70-200LB	Universal w/ Foot Drop Support	13"	6 Ea
M70-200LSB	Universal w. Foot Drop Support and Strap	13"	6 Ea
M70-200B	Bariatric w/Foot Drop Support and Strap	15"	4 Ea

Accessories

Product #	Description	Quantity
M70-100FW	Foam Wedge Only	6 Each