

Non-Adherent Super Absorbent Wound Dressings

- Innovative five-layer design
- Rounded and sealed edges eliminate fraying, loose threads, and linting
- Supports wound site protection from further trauma
- Soft pliability allows for range of motion exercises during recovery
- Able to apply saline or topical agents prior to use



Sofisorb® Specialty Shaped Dressings

- Tapered limb configurations conform to curved and hard-to-dress areas without unnecessary bulk
- Allows for customization for each patient's needs
- Available in over 50 sizes and configurations, including various anatomical shapes and pediatric sizes
- Supports ease of use

Five-Layer Construction

- 1 Non-adherent wound contact layer
- 2 Stay Dry™ Liner that permits passage of wound drainage to help keep periwound skin dry and prevent maceration at wound site
- 3 Center layer that absorbs wound drainage
- 4 Cellulose layer that wicks drainage horizontally along pad to increase dressing absorption capacity
- 5 Air-permeable backing for strength and durability



SOFSORB® Dressing 5-Layer Design

Sofsb® Specialty Absorbent Dressings

Tapered Limb

Item #	Description	Packaging	Quantity
46-111	XS (15" length)	1/Pk, 20 Pk/Cs	20/Cs
46-112	S (18" length)	1/Pk, 20 Pk/Cs	20/Cs
46-113	M (21" length)	1/Pk, 20 Pk/Cs	20/Cs
46-114	ML (23" length)	1/Pk, 15 Pk/Cs	15/Cs
46-115	L (26" length)	1/Pk, 15 Pk/Cs	15/Cs
46-116	XL (29" length)	1/Pk, 15 Pk/Cs	15/Cs



Standard Sizes

Item #	Description	Packaging	Quantity
46-109	3" x 3" (Round)	1/Pk, 30 Pk/Bx, 2 Bx/Cs	60/Cs
46-165	3" x 3" (Square)	1/Pk, 30 Pk/Bx, 2 Bx/Cs	60/Cs
46-162	4" x 4"	1/Pk, 50/Cs	50/Cs
46-101	4" x 6"	1/Pk, 30 Pk/Bx, 4 Bx/Cs	120/Cs
46-132	4" x 6" w/slit for trach/tube	1/Pk, 25 Pk/Bx, 4 Bx/Cs	100/Cs
46-106	4" x 9"	1/Pk, 25 Pk/Bx, 4 Bx/Cs	100/Cs
46-102	6" x 9"	1/Pk, 20 Pk/Bx, 4 Bx/Cs	80/Cs
46-103	9" x 15"	1/Pk, 30 Pk/Cs	30/Cs
46-104	15" x 18"	1/Pk, 20 Pk/Cs	20/Cs
46-105	15" x 24"	1/Pk, 15 Pk/Cs	15/Cs
46-166	20" x 28"	1/Pk, 20 Pk/Cs	20/Cs



Sofsb® Specialized Shape Absorbent Dressings

Item #	Description	Packaging	Quantity
46-121	Arm Dressing w/ties (Adult)	1/Pk, 20 Pk/Cs	20/Cs
46-153	Arm/Shoulder Dressing w/ties (Adult)	1/Pk, 20 Pk/Cs	20/Cs
46-127	Elbow Dressing	1/Pk, 20 Pk/Cs	20/Cs
46-122	Leg Dressing w/ties (Adult)	1/Pk, 20 Pk/Cs	20/Cs
46-155	Leg/Hip Dressing w/ties (Adult)	1/Pk, 10 Pk/Cs	10/Cs
46-133	Pediatric Chest/Vest (Small)	1/Pk, 20 Pk/Cs	20/Cs
46-134	Chest/Vest (Medium)	1/Pk, 20 Pk/Cs	20/Cs
46-123	Chest/Vest (Large)	1/Pk, 20 Pk/Cs	20/Cs
46-152	Chest/Vest (Extra-large)	1/Pk, 20 Pk/Cs	20/Cs
46-158	Jacket (Medium)	1/Pk, 20 Pk/Cs	20/Cs
46-151	Jacket (Large)	1/Pk, 10 Pk/Cs	10/Cs
46-128	Scalp Dressing	1/Pk, 20 Pk/Cs	20/Cs
46-129	Face Dressing (Adult)	1/Pk, 20 Pk/Cs	20/Cs
46-138	Pediatric Buttocks (Small)	1/Pk, 20 Pk/Cs	20/Cs
46-137	Buttocks (Medium)	1/Pk, 20 Pk/Cs	20/Cs
46-130	Buttocks (Large)	1/Pk, 15 Pk/Cs	15/Cs
46-159	Neck (Adult)	1/Pk, 20 Pk/Cs	20/Cs
46-156	Pediatric Hand/Glove (Small)	1/Pk, 20 Pk/Cs	20/Cs
46-131	Hand/Glove (Medium)	1/Pk, 20 Pk/Cs	20/Cs
46-157	Hand/Glove (Large)	1/Pk, 20 Pk/Cs	20/Cs
46-124	Pediatric Boot/Foot (Small)	1/Pk, 20 Pk/Cs	20/Cs
46-125	Boot/Foot (Medium)	1/Pk, 20 Pk/Cs	20/Cs
46-126	Boot/Foot (Large)	1/Pk, 20 Pk/Cs	20/Cs
46-160	Boot/Foot (Extra-large)	1/Pk, 10 Pk/Cs	10/Cs



**To learn more about our products or to place an order,
Call 1-888-938-7828 or visit deroyal.com.**