

## Features & Benefits:

- Polystyrene bead fill restrains flexion of the arm to help prevent self-extubation, IV dislodgement, and/or pulling at tubes, drains, and lines
- Non-metal construction allows for use in MRI or X-ray
- 3-strap hook and loop configuration for quick and easy application
- Wipeable vinyl exterior for easy cleaning
- Soft, breathable mesh interior for patient comfort
- Available in pediatric and adult sizes



Product No.	Size	Dimensions	UOM
M8035-I	Infant	3" L x 4" W	1 Ea
M8035-P	Pediatric	7" L x 7" W	1 Ea
M8035-C	Child	10" L x 12" W	1 Ea
M8035-S	Small	11" L x 14" W	1 Ea
M8035-M	Medium	13" L x 18" W	1 Ea
M8035-L	Large	14" L x 25" W	1 Ea