



Wound Care

SOLUTIONS

Table of Contents

Burn & General Wound Care

EXUDATE MANAGEMENT

Super Absorbents

- Sofsorb® Non-Adherent Absorbent Dressings7
- Sofsorb® Non-Adherent Specialty Shape Absorbent Dressings.....8
- Samson™ Non-Adherent Permeable / Impermeable Sheets.....12

Specialized Absorbents

- Specialized Shape Non-Adherent Gauze Dressings.....13
- Fluftex® High Performance Gauze Dressings.....15
- One-Plus Gauze Dressings.....16

Standard Absorbents

- Fine Mesh Gauze17
- Wide Mesh Gauze18
- Muslin Gauze / Gauze Fluff Pads20
- Cotton Rolls / Pads / Combine Dressings.....21

Non-Adherent Absorbents

- Multipad™ Dressings.....22

SECUREMENT / CLOSURE

- Stretch Net™ Elastic Retention Tubular Dressings23
- Apex™ Conforming Bandage.....24
- Bias Stockinette Absorbent Outer Wraps25
- Episeal® Wound Closure Strips.....26

COMPRESSION

- DeWrap™ 3-Layer Compression System27
- Unna Boot Zinc Oxide Impregnated Gauze28

JET LAVAGE DEBRIDEMENT

- JetOx™ Debridement System29

Advanced Wound Care

ANTIMICROBIAL DRESSINGS

Algidex Ag®+ Dressings31

- Algidex Ag®+ Foam Dressings.....32
- Algidex Ag®+ Pre-Cut Trach.....32
- Algidex Ag®+ Hydrogel / Packing Gauze / Wet Gauze Dressings33
- Algidex Ag®+ I.V. Dressings34

Covaderm® Ag+ Composite Dressings

- Covaderm® Ag+ Bordered Dressings35
- Covaderm® Ag+ Vascular Access Dressings (VAD)35

SPECIALTY WOUND TECHNOLOGY DRESSINGS

Multidex® Maltodextrin Wound Dressings

- Multidex® Gel.....36
- Multidex® Powder36

Table of Contents (cont.)

Advanced Wound Care (cont.)

FOAMS – ABSORBENT MULTI-LAYER

- Polyderm™ GTL⁵ 5-Layer Silicone Foam Dressings37
- Polyderm™ GTLite Silicone Foam Dressings37
- Polyderm™ Plus High Absorbency Foam Dressings41

COMPOSITE / ISLAND DRESSINGS

- Covaderm® Adhesive Wound Dressings42
- Covaderm® Thin Dressings43
- Covaderm Plus® Bordered / Vascular Access44

CONTACT LAYERS

- Dermanet® Non-Adherent Contact Layer Dressings.....45
- Dermanet® GTL Silicone Contact Layer Dressings.....46
- Dermanet® Ag+ / Ag+ Border / Ag+ Plus Border Absorbent Dressings.....47

FILM DRESSINGS

- Transeal® Transparent Film Dressings48
- Transeal® Plus Film Dressings49
- Transeal® TR / MAXX Film Dressings50

Negative Pressure Wound Therapy (NPWT)

INVENTORY MANAGEMENT

- Continuum® NPWT System Overview51
- Variable Pressure Therapy (VPT®) Technology53

TRADITIONAL NPWT

- Prospera® PRO-II™ NPWT System54
- Prospera® PRO NPWT Kits / Specialty Dressings.....55
- Prospera® Dressing Components / Accessories58

DISPOSABLE NPWT

- Prospera Flex® NPWT System60

NPWT COMPATIBLE DRESSINGS & ACCESSORIES

- Lumitech® Incisional Dressings.....62
- Incision Management Convenience Kits63

Appendix

- Healthcare Common Procedure Coding System64
- Trademark Attributions66
- References66



Burn and General Wound Care

Super Absorbents

Super Absorbent dressings are multi-layered dressings designed to manage highly exuding wounds. DeRoyal's Sofsorb® technology products and Samson™ sheets are available in a variety of shapes and sizes, allowing for customization to each patient's needs and supporting ease of use.

Sofsorb® Absorptive and Specialty Shaped Dressings

Sofsorb® is a highly absorbent, multi-layered 1-piece dressing. The dressing can be used dry or wet and secured in place by Stretch Net™ elastic bandage. Available in over 50 sizes and configurations, including various anatomical shapes and pediatric sizes.

Non-Adherent Super Absorbent Wound Dressings

- Innovative 5-layer design
- Rounded and sealed edges eliminate frayed edges, loose threads, and linting
- Tapered limb configurations conform to curved and hard-to-dress areas without unnecessary bulk
- Supports wound site protection from further trauma
- Soft pliability allows for range of motion exercises during recovery
- Appropriate for most wounds of various etiologies
- Can be used up to 7 days under compression for low to highly exudative environments
- Cost effective solution compared to traditional dressings which require frequent dressing changes
- Saline or topical agents may be applied prior to use



Sofsorb® 5-Layer Design

5 layer construction:

- 1 Non-adherent wound contact layer
- 2 StayDry™ Liner that permits passage of wound drainage to help keep periwound skin dry and prevent maceration at wound site
- 3 Center layer that absorbs wound drainage
- 4 Cellulose layer that wicks drainage horizontally along pad to increase dressing absorption capacity
- 5 Air-permeable backing for strength and durability

Exudate Management

Super Absorbents

Sofisorb® Absorbent Dressings

Standard Sizes

| Item # | Description | Packaging | Quantity |
|--------|----------------------------------|-------------------------|----------|
| 46-109 | 3" x 3" (Round) | 1/Pk, 30 Pk/Bx, 2 Bx/Cs | 60/Cs |
| 46-165 | 3" x 3" (Square) | 1/Pk, 30 Pk/Bx, 2 Bx/Cs | 60/Cs |
| 46-162 | 4" x 4" | 1/Pk, 50/Cs | 50/Cs |
| 46-101 | 4" x 6" | 1/Pk, 30 Pk/Bx, 4 Bx/Cs | 120/Cs |
| 46-132 | 4" x 6" w/ Slit for Trach / Tube | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |
| 46-106 | 4" x 9" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |
| 46-102 | 6" x 9" | 1/Pk, 20 Pk/Bx, 4 Bx/Cs | 80/Cs |
| 46-103 | 9" x 15" | 1/Pk, 30 Pk/Cs | 30/Cs |
| 46-104 | 15" x 18" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-105 | 15" x 24" | 1/Pk, 15 Pk/Cs | 15/Cs |
| 46-166 | 20" x 28" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-107 | 24" x 36" | 1/Pk, 20 Pk/Cs | 20/Cs |



Sofisorb® Specialty Shape Absorbent Dressings

Tapered Limb

| Item # | Description | Packaging | Quantity |
|--------|---------------------------|----------------|----------|
| 46-111 | Extra Small (15" Length) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-112 | Small (18" Length) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-113 | Medium (21" Length) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-114 | Medium Large (23" Length) | 1/Pk, 15 Pk/Cs | 15/Cs |
| 46-115 | Large (26" Length) | 1/Pk, 15 Pk/Cs | 15/Cs |
| 46-116 | Extra Large (29" Length) | 1/Pk, 15 Pk/Cs | 15/Cs |



Soforb® Specialized Shape Absorbent Dressings

Arm Dressings

| Item # | Description | Packaging | Quantity |
|--------|---|----------------|----------|
| 46-121 | Arm Dressing w/ Ties (Adult) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-153 | Arm / Shoulder Dressing w/ Ties (Adult) | 1/Pk, 20 Pk/Cs | 20/Cs |



Elbow Dressings

| Item # | Description | Packaging | Quantity |
|--------|----------------|----------------|----------|
| 46-127 | Elbow Dressing | 1/Pk, 20 Pk/Cs | 20/Cs |



Leg Dressings

| Item # | Description | Packaging | Quantity |
|--------|------------------------------------|----------------|----------|
| 46-122 | Leg Dressing w/ Ties (Adult) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-155 | Leg / Hip Dressing w/ Ties (Adult) | 1/Pk, 10 Pk/Cs | 10/Cs |



Exudate Management

Super Absorbents

Sofsb® Specialized Shape Absorbent Dressings

Chest Dressings / Vests

| Item # | Description | Packaging | Quantity |
|--------|--------------------------------|----------------|----------|
| 46-133 | Pediatric Chest / Vest (Small) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-134 | Chest / Vest (Medium) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-123 | Chest / Vest (Large) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-152 | Chest / Vest (Extra-large) | 1/Pk, 20 Pk/Cs | 20/Cs |



Jackets

| Item # | Description | Packaging | Quantity |
|--------|-----------------|----------------|----------|
| 46-158 | Jacket (Medium) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-151 | Jacket (Large) | 1/Pk, 10 Pk/Cs | 10/Cs |



Head / Scalp Dressings

| Item # | Description | Packaging | Quantity |
|--------|----------------|----------------|----------|
| 46-128 | Scalp Dressing | 1/Pk, 20 Pk/Cs | 20/Cs |



Face Dressings

| Item # | Description | Packaging | Quantity |
|--------|-----------------------|----------------|----------|
| 46-129 | Face Dressing (Adult) | 1/Pk, 20 Pk/Cs | 20/Cs |

Sofsb® Specialized Shape Absorbent Dressings

Buttocks Dressings

| Item # | Description | Packaging | Quantity |
|--------|----------------------------|----------------|----------|
| 46-138 | Pediatric Buttocks (Small) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-137 | Buttocks (Medium) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-130 | Buttocks (Large) | 1/Pk, 15 Pk/Cs | 15/Cs |



Neck Dressings

| Item # | Description | Packaging | Quantity |
|--------|--------------|----------------|----------|
| 46-159 | Neck (Adult) | 1/Pk, 20 Pk/Cs | 20/Cs |



Hand Dressings

| Item # | Description | Packaging | Quantity |
|--------|--------------------------------|----------------|----------|
| 46-156 | Pediatric Hand / Glove (Small) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-131 | Hand / Glove (Medium) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-157 | Hand / Glove (Large) | 1/Pk, 20 Pk/Cs | 20/Cs |



Foot Dressings

| Item # | Description | Packaging | Quantity |
|--------|-------------------------------|----------------|----------|
| 46-124 | Pediatric Boot / Foot (Small) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-125 | Boot / Foot (Medium) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-126 | Boot / Foot (Large) | 1/Pk, 20 Pk/Cs | 20/Cs |
| 46-160 | Boot / Foot (Extra-large) | 1/Pk, 10 Pk/Cs | 10/Cs |




Exudate Management

Super Absorbents

Samson™ Sheets Non-Adherent Permeable / Impermeable Sheets

DeRoyal offers a unique multi-layered burn sheet designed with several features:

- First layer is composed of a non-adherent contact layer to limit adherence to the wound area
- Second layer is composed of several sheets of compressed cellulose for exudate absorbency
- Third layer provides reinforcement to support use as a transfer or draw sheet
- Optional permeable or impermeable backing



Impermeable

| Item # | Description | Packaging | Quantity |
|---------|--------------------------|-----------------|----------|
| 10-6436 | 24" x 36" | 1/Pk, 24/Cs | 24/Cs |
| 10-8126 | 30" x 36" | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-8130 | 32" x 45" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-8134 | 32" x 72" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-8131 | 32" x 90" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-8135 | 16" x 22" w/ Cheesecloth | 1/Pk, 30 Pk/Cs | 30/Cs |
| 10-2824 | 24" x 36" w/ Cheesecloth | 1/Pk, 15 Pk/Cs | 15/Cs |
| 10-6536 | 24" x 36" w/ Cheesecloth | 1/Pk, 24 Pk/Cs | 24/Cs |

Permeable

| Item # | Description | Packaging | Quantity |
|---------|-------------|----------------|----------|
| 10-8133 | 32" x 72" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-8138 | 32" x 45" | 1/Pk, 20 Pk/Cs | 20/Cs |

Specialized Absorbents

Special shaped gauze dressings have the same function and benefits as regular gauze dressings but have been anatomically cut to fit a particular part of the body. These dressings are easy to apply, save nursing time, and increase patient comfort.

Specialized Shape Non-Adherent Gauze Dressings

Vests

| Item # | Description | Packaging | Quantity |
|-----------|--------------------------|----------------|----------|
| 10-9134 | 9" x 12" 15-ply (Infant) | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-9133 | 14" x 18" 15-ply (Child) | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6314 | 14" x 18" 24-ply (Child) | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6313 | 22" x 26" 24-ply (Teen) | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-9131 | 28" x 32" 15-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6447 | 30" x 32" 15-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-9138 | 30" x 32" 32-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-5514 | 30" x 36" 28-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-9137 | 30" x 52" 16-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-9137-L | 30" x 64" 16-ply | 1/Pk, 20 Pk/Cs | 20/Cs |



Arm Dressings

| Item # | Description | Packaging | Quantity |
|---------|------------------|----------------|----------|
| 10-1022 | 17" x 35" 16-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6437 | 19" x 15" x 28" | 1/Pk, 6 Pk/Cs | 6/Cs |



Exudate Management

Specialized Absorbents

Specialized Shape Non-Adherent Gauze Dressings

Leg Dressings

| Item # | Description | Packaging | Quantity |
|---------|------------------|----------------|----------|
| 10-6321 | 20" x 27" 24-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6306 | 36" x 36" 24-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-6438 | 29" x 24" x 40" | 1/Pk, 6 Pk/Cs | 6/Cs |



Face Dressings

| Item # | Description | Packaging | Quantity |
|---------|---|----------------|----------|
| 10-8173 | 12" x 14" 10-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6467 | 16" x 24" 20-ply w/ Neck Flaps & Tube Slit | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6468 | 12" x 18" 20-ply w/ Neck Flaps & Tube Slit (Pediatric) | 1/Pk, 25 Pk/Cs | 25/Cs |



Buttocks Dressings

| Item # | Description | Packaging | Quantity |
|---------|------------------|----------------|----------|
| 10-5515 | 22" x 60" 28-ply | 1/Pk, 25 Pk/Cs | 25/Cs |



Fluftex® Gauze Dressings

Fluftex® products are a line of gauze dressings that have been pre-washed, fluff-dried, and processed for superior absorption and softness.

- Specially finished or folded edges to reduce lint
- Ideal for wound management, preparation, packing, cleansing and bandaging
- Available in sponges, rolls, or large pad options



Sponges

| Item # | Description | Packaging | Quantity |
|--------|---|--------------------------|----------|
| 11-001 | 6" x 6¾", Medium Sterile 2's in Soft Pouch | 2/Pk, 20 Pk/Bx, 12 Bx/Cs | 480/Cs |
| 11-003 | 6" x 6¾", Medium Sterile 5's in Soft Pouch | 5/Pk, 10 Pk/Bx, 12 Bx/Cs | 600/Cs |
| 11-004 | 6" x 6¾", Medium Sterile 10's in Plastic Tray | 10/Pk, 48 Pk/Cs | 480/Cs |
| 11-005 | 6" x 6¾", Medium Non-sterile, Bulk | 100/Pk, 6 Pk/Cs | 600/Cs |
| 11-006 | 7¾" x 8¾", Large Non-sterile, Bulk | 100/Pk, 4 Pk/Cs | 400/Cs |

Rolls

| Item # | Description | Packaging | Quantity |
|--------|-------------------------------------|-----------------|----------|
| 11-009 | 2¼" x 3 yds, Sterile in Soft Pouch | 1/Pk, 96 Pk/Cs | 96/Cs |
| 11-008 | 3⅝" x 3⅝ yds, Sterile in Soft Pouch | 1/Pk, 96 Pk/Cs | 96/Cs |
| 11-012 | 3⅝" x 3⅝ yds, Non-sterile, Bulk | 96/Cs | 96/Cs |
| 11-007 | 4½" x 4⅓ yds, Sterile in Soft Pouch | 1/Pk, 100 Pk/Cs | 100/Cs |
| 11-011 | 4½" x 4⅓ yds, Non-sterile, Bulk | 100/Cs | 100/Cs |

Pads

| Item # | Description | Packaging | Quantity |
|---------|-----------------|----------------|----------|
| 32-248 | 9" x 12" 15-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-5139 | 9" x 28" 15-ply | 2/Pk, 20 Pk/Cs | 40/Cs |

Exudate Management

Specialized Absorbents

One-Plus Gauze Dressings

A One-Plus dressing is made by serging 1 layer of Fine Mesh to several layers of Wide Mesh.

- 2-sided design for Wide Mesh or Fine Mesh gauze wound contact options
- Wide Mesh side supports debridement needs
- Fine Mesh side limits wound adherence and supports exudate management
- Eliminates the excessive time required to pre-cut and layer gauze
- Available in a variety of sizes and ply options



Gauze Dressings

| Item # Quantity | Description | Packaging | |
|--------------------|------------------|----------------|--------|
| 10-5500 | 18" x 18" 6-ply | 2/Pk, 50 Pk/Cs | 100/Cs |
| 10-5510 | 36" x 36" 17-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 15-018 | 36" x 36" 18-ply | 1/Pk, 10 Pk/Cs | 10/Cs |
| 15-025 | 18" x 18" 6-ply | 4/Pk, 25 Pk/Cs | 100/Cs |

Fine Mesh Gauze Dressings

Fine Mesh gauze is less adherent to wound surfaces and can be used to cover donor sites and skin grafts.

- Available as sheets and rolls
- Topical solution or antibacterial cream may be applied prior to application



Sheets

| Item # | Description | Packaging | Quantity |
|---------|-------------------------------|-----------------|----------|
| 10-6430 | 2" x 3" 1-ply | 5/Pk, 50 Pk/Cs | 250/Cs |
| 10-038 | 3" x 8" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-6393 | 3" x 36" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-039 | 4" x 8" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-410 | 4" x 10" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-4110 | 4" x 10" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-425 | 4" x 25" 1-ply | 2/Pk, 100 Pk/Cs | 200/Cs |
| 10-912 | 9" x 12" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-918 | 9" x 18" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-1318 | 13" x 18" 1-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-6367 | 18" x 18" 1-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-1804 | 18" x 18" 1-ply | 4/Pk, 20 Pk/Cs | 80/Cs |
| 10-1818 | 18" x 18" 1-ply (Non-sterile) | 4/Pk, 20 Pk/Cs | 80/Cs |
| 10-4124 | 24" x 24" 12-ply | 1/Pk, 30 Pk/Cs | 30/Cs |

Rolls

| Item # | Description | Packaging | Quantity |
|---------|-------------------|----------------|----------|
| 10-9130 | 2" x 5 yds 1-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 32-033 | 3" x 5 yds 1-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-4360 | 4" x 10 yds 1-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-6360 | 6" x 10 yds 1-ply | 1/Pk, 50 Pk/Cs | 50/Cs |

Exudate Management

Standard Absorbents

Wide Mesh Unfinished Edge Gauze Dressings

Wide Mesh gauze dressings are used for absorption, padding and debridement.

- Unique die-cutting limits lint and loose threading
- Available in a variety of sizes and ply options
- Topical solution or antibacterial cream may be applied prior to application



Unfinished Edge Gauze Dressings

| Item # | Description | Packaging | Quantity |
|---------|-------------------------------|-----------------|----------|
| 10-5121 | 4" x 72" 12-ply (Non-sterile) | 10/Pk, 5 Pk/Cs | 50/Cs |
| 10-4113 | 9" x 36" 24-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-6401 | 9" x 36" 4-ply | 12/Pk, 50 Pk/Cs | 600/Cs |
| 32-021 | 9" x 72" 4-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5126 | 9" x 72" 12-ply (Non-sterile) | 10/Pk, 5 Pk/Cs | 50/Cs |
| 10-6440 | 9" x 72" 4-ply | 2/Pk, 25 Pk/Cs | 50/Cs |
| 10-6441 | 9" x 72" 4-ply | 5/Pk, 20 Pk/Cs | 100/Cs |
| 10-767 | 12" x 12" 12-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-1001 | 18" x 18" 5-ply | 6/Pk, 20 Pk/Cs | 120/Cs |
| 10-5173 | 18" x 18" 6-ply | 2/Pk, 50 Pk/Cs | 100/Cs |
| 10-5111 | 18" x 18" 8-ply | 4/Pk, 50 Pk/Cs | 200/Cs |
| 10-8171 | 18" x 18" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5107 | 18" x 18" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5167 | 18" x 18" 12-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-6400 | 18" x 18" 16-ply | 6/Pk, 50 Pk/Cs | 300/Cs |
| 10-1820 | 18" x 18" 20-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5106 | 18" x 18" 50-ply | 1/Pk, 36 Pk/Cs | 36/Cs |
| 10-6537 | 18" x 18" 100-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-5161 | 18" x 36" 4-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5109 | 18" x 36" 4-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-5105 | 18" x 36" 8-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-766 | 18" x 36" 12-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6566 | 18" x 36" 14-ply | 1/Pk, 25 Pk/Cs | 25/Cs |

Wide Mesh Unfinished Edge Gauze Dressings

| Item # | Description | Packaging | Quantity |
|---------|------------------|-----------------|----------|
| 10-5169 | 18" x 36" 20-ply | 1/Pk, 30 Pk/Cs | 30/Cs |
| 10-8172 | 18" x 60" 8-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-5164 | 24" x 30" 20-ply | 1/Pk, 30 Pk/Cs | 30/Cs |
| 10-1500 | 24" x 36" 8-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-6579 | 24" x 36" 8-ply | 4/Pk, 15/Cs | 60/Cs |
| 10-5130 | 36" x 36" 1-ply | 2/Pk, 50 Pk/Cs | 100/Cs |
| 10-6453 | 36" x 36" 1-ply | 3/Pk, 100 Pk/Cs | 300/Cs |
| 10-6108 | 36" x 36" 10-ply | 1/Pk, 30 Pk/Cs | 30/Cs |
| 10-768 | 36" x 36" 12-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 15-003 | 12" x 12" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 15-009 | 18" x 18" 10-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 15-011 | 18" x 18" 50-ply | 1/Pk, 30 Pk/Cs | 30/Cs |
| 15-019 | 36" x 36" 18-ply | 1/Pk, 20 Pk/Cs | 20/Cs |
| 15-023 | 36" x 36" 10-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-6458 | 18" x 184" 2-ply | 1/Pk, 25 Pk/Cs | 25/Cs |

Wide Mesh Finished Edge Gauze Dressings

DeRoyal offers a variety of wide mesh gauze dressings with finished edges.

- Available in several sizes, shapes and ply options
- Available in serged or sewn edges



Finished Edge Gauze Dressings

| Item # | Description | Packaging | Quantity |
|---------|--|-----------------|----------|
| 10-9000 | 18" x 18" 5-ply | 6/Pk, 20 Pk/Cs | 120/Cs |
| 10-8143 | 18" x 18" 10-ply | 3/Pk, 20 Pk/Cs | 60/Cs |
| 10-6322 | 18" x 18" 12-ply (Dermanet® Contact Layer) | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-8144 | 18" x 18" 15-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-932 | 18" x 18" 18-ply | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-5181 | 18" x 36" 32-ply | 1/Pk, 25 Pk/Cs | 25/Cs |
| 10-6300 | 24" x 36" 16-ply | 1/Pk, 30 Pk/Cs | 30/Cs |
| 10-6448 | 36" x 28" 15-ply | 1/Pk, 20 Pk/Cs | 20/Cs |

Note: Finished edges cannot be separated

Exudate Management

Standard Absorbents

Muslin Gauze Absorbent Towel

Muslin Gauze is an absorbent towel comprised of 3 layers of material: 1-ply cheesecloth, 4-ply wide mesh gauze, and 3-ply fine mesh gauze with serged edges.



Muslin Gauze

| Item # | Description | Packaging | Quantity |
|---------|-------------|----------------|----------|
| 10-6388 | 18" x 27" | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-6548 | 18" x 35" | 1/Pk, 50 Pk/Cs | 50/Cs |

Gauze Fluff Absorbent Pads

Gauze fluff pads provides absorptive padding that is appropriate for:

- Heavily exudating wounds (e.g. burn)
- Wiping and debriding wound areas



Gauze Fluff

| Item # | Description | Packaging | Quantity |
|---------|-------------------------------|-----------------|----------|
| 31-466 | 15" x 15" 2-ply | 15/Pk, 20 Pk/Cs | 300/Cs |
| 10-8153 | 18" x 12" 1-ply (Non-sterile) | 500/Cs | 500/Cs |
| 10-8156 | 18" x 36" 1-ply | 4/Pk, 20 Pk/Cs | 80/Cs |
| 10-8154 | 18" x 36" 1-ply (Non-sterile) | 500/Cs | 500/Cs |
| 10-8155 | 36" x 36" 1-ply (Non-sterile) | 500/Cs | 500/Cs |
| 10-8150 | 36" x 36" 1-ply | 1/Pk, 250 Pk/Cs | 250/Cs |
| 10-8157 | 36" x 36" 1-ply | 4/Pk, 20 Pk/Cs | 80/Cs |
| 31-980 | 36" x 36" 1-ply | 5/Pk, 50 Pk/Cs | 250/Cs |
| 10-8152 | 36" x 36" 2-ply | 1/Pk, 100 Pk/Cs | 100/Cs |
| 10-8158 | 36" x 36" 2-ply | 4/Pk, 20 Pk/Cs | 80/Cs |

Cotton Rolls / Pads / Combine Dressings

Cotton rolls, cotton pads, and combine dressings offer superior absorption capabilities and light compression. The sizes can be customized to fit individual needs.



Rolls

| Item # | Description | Packaging | Quantity |
|---------|---------------------|----------------|----------|
| 9866-01 | ½ lb Roll 6" x 13' | 1/Pk, 10 Pk/Cs | 10/Cs |
| 9866-00 | 1 lb Roll 12" x 11' | 1/Pk, 10 Pk/Cs | 10/Cs |

Pads

| Item # | Description | Packaging | Quantity |
|--------|-------------|-----------------|----------|
| 31-216 | 3" x 3" | 1/Pk, 100 Pk/Cs | 100/Cs |
| 31-091 | 4" x 4" | 1/Pk, 100 Pk/Cs | 100/Cs |
| 31-783 | 6" x 6" | 1/Pk, 50 Pk/Cs | 50/Cs |
| 31-475 | 6" x 8" | 1/Pk, 50 Pk/Cs | 50/Cs |
| 31-648 | 6" x 12" | 1/Pk, 25 Pk/Cs | 25/Cs |

Dressings

| Item # | Description | Packaging | Quantity |
|---------|-------------|----------------|----------|
| 10-5140 | 8" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 10-5141 | 8" x 10 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 31-407 | 8" x 24" | 1/Pk, 50 Pk/Cs | 50/Cs |
| 10-5143 | 8" x 30" | 1/Pk, 30 Pk/Cs | 30/Cs |

Exudate Management

Non-Adherent Absorbent Dressings

Non-Adherent Absorbent Dressings

Multipad™ dressings are non-adherent, highly absorbent wound dressings. The design is composed of a non-woven rayon pad between 2 non-adhering polypropylene contact layers.

Multipad™ Non-Adherent Wound Dressings

- A multi-layered, highly absorbent wound dressing
- Helps protect fragile granulation tissue and prevents wound trauma
- Non-linting
- Use in place of gauze and tape when more absorbency is needed
- Can be cut for oddly shaped wounds



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-011 | 2" x 2" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs | A6251 |
| 46-012 | 4" x 4" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs | A6251 |
| 46-013 | 4" x 8" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6252 |
| 46-014 | 7½" x 7½" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6253 |

Securement / Closure

DeRoyal offers a variety of dressing securement and wound closure options. Securement options include DeRoyal's proprietary Stretch Net™, Apex™, and Bias Stockinette dressing products. The Episeal™ wound closure strips provide clinicians with wound and incision closure options.

Stretch Net™ Tubular Elastic Retention Dressings

- Tubular elastic retention bandage that holds dressings in place without adhesives or tapes
- Allows easy wound inspection and quick application for wounds with frequent dressing changes
- Provides maximum aeration to wound
- Ideal for patients with fragile skin or adhesive allergies



Stretch Net™ Product Options (10 Yards / 25 Yards)

| Standard | | Thin | | Stretch Net™ Size | Approximate Circumference | Suggested Application | HCPCS |
|--------------------|--------------------|-------------------------|----|-------------------|---------------------------|------------------------|-------|
| 10 Yards Product # | 25 Yards Product # | 25 Yards Thin Product # | | | | | |
| 10-7401 | 10-7421 | 15-7001 | 1 | | 3" to 8" | Fingers, Toes | A6457 |
| 10-7402 | 10-7422 | 15-7002 | 2 | | 8" to 11" | Fingers, Wrist | A6457 |
| 10-7403 | 10-7423 | 15-7003 | 3 | | 11" to 16" | Hand, Elbow, Foot | A6457 |
| 10-7404 | 10-7424 | 15-7004 | 4 | | 16" to 22" | Hand, Elbow | A6457 |
| 10-7405 | 10-7425 | 15-7005 | 5 | | 22" to 26" | Calf, Thigh, Head | A6457 |
| 10-7406 | 10-7426 | 15-7006 | 5½ | | 26" to 35" | Small Chest, Axilla | A6457 |
| 10-7407 | 10-7427 | 15-7007 | 6 | | 35" to 42" | Abdomen, Thigh | A6457 |
| 10-7408 | 10-7428 | 15-7008 | 7 | | 42" to 52" | Chest, Axilla, Abdomen | A6457 |
| 10-7409 | 10-7429 | 15-7009 | 8 | | 52" to 58" | Chest, Axilla, Abdomen | A6457 |
| 10-7410 | 10-7430 | 15-7010 | 9 | | 58" to 64" | Chest, Axilla, Abdomen | A6457 |
| 10-7411 | 10-7431 | 15-7011 | 10 | | 64" to 72" | Chest, Axilla, Abdomen | A6457 |
| 10-7412 | 10-7432 | 15-7012 | 11 | | 72" to 80" | Extra Large Torso | A6457 |
| 10-7413 | 10-7433 | 15-7013 | 12 | | 80" to 86" | Extra Large Torso | A6457 |

Apex™ Conforming Bandage

Apex™ conforming bandages are low linting, 1-ply, rayon/polyester blend dressings.

- Conforms easily to body contours
- Self-Adherent
- Maintains conformability when wet
- Absorbent
- Soft and gentle to the skin
- Limits excess bulk



Sterile

| Product # | Description | Packaging | Quantity | HCPs |
|-----------|-------------|-----------------|----------|-------|
| 11-6991 | 1" x 60" | 1/Pk, 144 Pk/Cs | 144/Cs | A6445 |
| 11-6992 | 2" x 65" | 1/Pk, 96 Pk/Cs | 96/Cs | A6445 |
| 11-6993 | 3" x 75" | 1/Pk, 72 Pk/Cs | 72/Cs | A6446 |
| 11-6994 | 4" x 75" | 1/Pk, 72 Pk/Cs | 72/Cs | A6446 |
| 11-6996 | 6" x 85" | 1/Pk, 48 Pk/Cs | 48/Cs | A6447 |

Non-sterile

| Product # | Description | Packaging | Quantity | HCPs |
|-----------|-------------|----------------|----------|-------|
| 11-6981 | 1" x 60" | 24/Bg, 4 Bg/Cs | 96/Cs | A6445 |
| 11-6982 | 2" x 65" | 12/Bg, 8 Bg/Cs | 96/Cs | A6445 |
| 11-6983 | 3" x 75" | 12/Bg, 8 Bg/Cs | 96/Cs | A6446 |
| 11-6984 | 4" x 75" | 12/Bg, 8 Bg/Cs | 96/Cs | A6446 |
| 11-6986 | 6" x 85" | 6/Bg, 8 Bg/Cs | 48/Cs | A6447 |

Bias Stockinette Absorbent Outer Wraps

- 100% cotton bias stockinette rolls
- Used as outer wrap over other dressings
- Absorbent and provides extra support and comfort to the patient



Absorbent Outer Wrap

| Item # | Description | Packaging | Quantity |
|---------|--|----------------|----------|
| 9800-24 | 2" x 4 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-25 | 2" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-33 | 3" x 3 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-34 | 3" x 4 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-35 | 3" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-43 | 4" x 3 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-44 | 4" x 4 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-45 | 4" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-46 | 4" x 6 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-63 | 6" x 3 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-64 | 6" x 4 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-65 | 6" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-66 | 6" x 6 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-84 | 8" x 4 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9800-85 | 8" x 5 yds | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9804-56 | 4" x 56" | 1/Pk, 20 Pk/Cs | 20/Cs |
| 9844-43 | 4" x 3 yds (Non-sterile) | 4/Pk, 25 Pk/Cs | 100/Cs |
| 9833-63 | 6" x 3 yds (High-density, Non-sterile) | 4/Pk, 25 Pk/Cs | 100/Cs |
| 9866-63 | 6" x 3 yds (High-density, Non-sterile) | 4/Pk, 25 Pk/Cs | 100/Cs |
| RB2 | 2" x 50 yds | 1 RI | 1/Cs |
| RB3 | 3" x 50 yds | 1 RI | 1/Cs |
| RB4 | 4" x 50 yds | 1 RI | 1/Cs |
| RB4-3 | 4" x 3 yds | 1 RI | 1/Cs |
| RB4-5 | 4" x 5 yds | 1 RI | 1/Cs |
| RB5-3 | 5" x 3 yds | 1 RI | 1/Cs |
| RB6 | 6" x 50 yds | 1 RI | 1/Cs |
| RB6-3 | 6" x 3 yds | 1 RI | 1/Cs |
| RB6-5 | 6" x 5 yds | 1 RI | 1/Cs |

Securement / Closure

Closure

Episeal® Wound Closure Strips

- Wound closure strips incorporate a non-irritating adhesive coating
- Provides secure attachment without the need for additional adhesive products (i.e. benzoin tincture)
- Rounded corners limit edge lift and supports securement
- High breathability and tensile strength
- Flexible and conformable



Wound Closure Strips

| Product # | Description | Strips/Pack | Packaging | Quantity | HCPCS |
|-----------|---------------|-------------|-------------------|-----------|-------|
| 46-211 | 1/8" x 3" | 5 | 50 Pk/Bx, 4 Bx/Cs | 200 Pk/Cs | A4452 |
| 46-212 | 1/4" x 3" | 3 | 50 Pk/Bx, 4 Bx/Cs | 200 Pk/Cs | A4452 |
| 46-213 | 1/4" x 4" | 10 | 50 Pk/Bx, 4 Bx/Cs | 200 Pk/Cs | A4452 |
| 46-214 | 1/2" x 4" | 6 | 50 Pk/Bx, 4 Bx/Cs | 200 Pk/Cs | A4452 |
| 46-215 | 1/4" x 1 1/2" | 6 | 50 Pk/Bx, 4 Bx/Cs | 200 Pk/Cs | A4452 |
| 46-216 | 1" x 5" | 4 | 50 Pk/Bx, 4 Bx/Cs | 100 Pk/Cs | A4452 |

Compression

The DeRoyal system offers specialized compression wraps designed to aid in lower extremity wound healing. The DeWrap™ system is a 3-layer graduated compression system. Unna Boot is a specialized, zinc oxide impregnated gauze wrap. Both the DeWrap™ system and the Unna Boot are comfortable, conformable, and easy to apply.

DeWrap™ 3-Layer Compression System

- 3-layer component compression system that delivers sustained, graduated compression for lower extremities
- Reduces venous hypertension and aids in the healing of venous leg ulcerations
- Application instructions available to demonstrate how to achieve different levels of compression
- Not made with natural rubber latex



3-Layer Compression System:

- 1** Absorptive layer is absorbent, conformable, and protective
- 2** Elastic layer provides breathability and graduated pressure
- 3** Cohesive layer is self-adhering to secure bandage system and provides graduated compression for up to 7-day wear time

DeWrap™ 3-Layer Compression System

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|----------------|---------------|----------|--|
| 46-333 | 3-Layer System | 1/Pk, 8 Pk/Cs | 8/Cs | Layer 1: A6441 Layer 2: A6451 Layer 3: A6454 |

Compression

Zinc Oxide Impregnated Gauze

Unna Boot Zinc Oxide Impregnated Gauze

- Gauze impregnated with non-hardening zinc oxide paste
- Easily removed with blunt bandage scissors
- Convenient pack available to include cohesive wraps
- Multiple sizes and configurations available



Available Layers:

- 1 Zinc Oxide Impregnated Gauze
- 2 Cohesive Layer (Optional)

Zinc Impregnated Gauze

| Product # | Description | Packaging | Quantity |
|-----------|--|----------------|----------|
| 46-334 | 3" x 10 yds | 1/Pk, 12 Pk/Cs | 12/Cs |
| 46-335 | 4" x 10 yds | 1/Pk, 12 Pk/Cs | 12/Cs |
| 46-336 | 3" x 10 yds Unna Boot 4" x 6½ yds Cohesive Wrap | 1/Pk, 12 Pk/Cs | 12/Cs |
| 46-337 | 4" x 10 yds Unna Boot 4" x 6½ yds Cohesive Wrap | 1/Pk, 12 Pk/Cs | 12/Cs |

Jet Lavage Debridement

JetOx™ Debridement System

Jet Lavage Debridement

JetOx™-ND initiates wound cleansing and debridement through compressed air and oxygen. The JetOx™-ND system administers a cool, gentle stream that creates a desensitizing effect.

JetOx™-ND Jet Lavage Wound Cleansing & Debridement System

- Simple, user-controlled pressure
- Quickly and efficiently cleans and debrides wounds
- Delivers a cool, gentle fluid jet stream
- Microdrops between 5-200 microns accelerated up to 200 meters per second
- For patients sensitive to other forms of debridement
- Compatible with standard gas outlets (portable and wall) and saline dispensers
- Ideal for bedside application
- Fully disposable, limiting cleanup of the work environment



JetOx™-ND Lavage Product Options

| Product # | Description | Packaging | Quantity |
|-----------|---|----------------|----------|
| 48-500 | JetOx™-ND Wand | 1/Pk, 25 Pk/Cs | 25/Cs |
| 48-601 | JetOx™-ND Fit all Tubing Sets | 1/Pk, 25 Pk/Cs | 25/Cs |
| 48-444 | JetOx™-ND Mini-pack w/ 5 Wands, 5 Fit all Tubing Sets | 1/Pk, 5 Pk/Cs | 5/Cs |



Advanced Wound Care Dressings

Antimicrobial Dressings

Algidex Ag®+



Algidex Ag®+ Antimicrobial Wound Dressings

Algidex Ag®+ is a unique combination of silver, alginate, and maltodextrin designed to support an optimal moist wound environment.

- Provides quick and sustained antimicrobial activity against a broad spectrum of microbial organisms for up to 7 days¹
- Does not require activation or the use of sterile distilled water
- Up to 7 day wear time limits frequent dressing changes
- Inhibits bacterial cell division which minimizes the development of resistance*
- Gentle and non-irritating

* Data on File

¹ Broad spectrum represented in testing by the following microbial organisms: Escherichia coli, Pseudomonas aeruginosa, Staphylococcus aureus, methicillin-resistant Staphylococcus aureus (MRSA), and Candida albicans

Antimicrobial Dressings

Algidex Ag[®]+

Algidex Ag[®]+ Foam Dressings

- Non-adherent antimicrobial foam dressing
- Used for low to moderate draining wounds
- Absorptive to help reduce risk for periwound maceration
- Thick, soft construction for conformability and comfort



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------------|-------------------------|----------|-------|
| 46-AG22 | 2" x 2" Foam Back | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6209 |
| 46-AG44 | 4" x 4" Foam Back | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6209 |
| 46-AG45 | 4" x 5" Foam Back | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6210 |
| 46-AG66 | 6" x 6" Foam Back | 1/Pk, 5 Pk/Bx, 5 Bx/Cs | 25/Cs | A6210 |
| 46-AG88 | 8" x 8" Foam Back | 1/Pk, 5 Pk/Bx, 5 Bx/Cs | 25/Cs | A6211 |

Algidex Ag[®]+ Pre-Cut Trach

- Pre-cut for easy application and removal
- Non-adherent antimicrobial foam dressing
- Thick, soft construction for conformability and comfort



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-----------------------------------|-------------------------|----------|-------|
| 46-TD22 | 2" x 2" Pre-cut Trach Dressing | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6209 |
| 46-TD44 | 4" x 4" Pre-cut Trach Dressing | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6209 |

Antimicrobial Dressings

Algidex Ag[®]+

Algidex Ag[®]+ Hydrogel

- Hydrogel ionic silver dressing
- Provides antimicrobial benefits for dry to low draining wounds
- Supports autolytic debridement
- For use on infected or non-infected wounds



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-G10 | 1 fl oz | 1/Pk, 3 Pk/Bx, 10 Bx/Cs | 30/Cs | A6248 |

Algidex Ag[®]+ Packing Gauze

- Ideal for deep or tunneling wounds
- Pre-moistened for ease of placement
- Specialty woven packing gauze resistant to unraveling, linting, or fraying



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|----------------|----------|-------|
| 46-PG12 | ½" x 5 yds | 1/Bt, 12 Bt/Cs | 12/Cs | A6407 |
| 46-PG14 | ¼" x 5 yds | 1/Bt, 12 Bt/Cs | 12/Cs | A6407 |

Algidex Ag[®]+ Wet Gauze

- Ideal for dry to minimally draining wounds
- Specialty woven gauze is resistant to unraveling, linting, or fraying
- Pre-moistened for ease of placement
- Provides quick and sustained antimicrobial activity against a broad spectrum of microbial organisms for up to 7 days¹



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-GZ22 | 2" x 2" | 1/Pk, 20 Pk/Bx, 5 Bx/Cs | 100/Cs | A6231 |
| 46-GZ44 | 4" x 4" | 1/Pk, 20 Pk/Bx, 5 Bx/Cs | 100/Cs | A6231 |

¹ Broad spectrum represented in testing by the following microbial organisms: Escherichia coli, Pseudomonas aeruginosa, Staphylococcus aureus, methicillin-resistant Staphylococcus aureus (MRSA), and Candida albicans

Antimicrobial Dressings

Algidex Ag[®]+

Algidex Ag[®]+ I.V. Dressings

- Provides quick and sustained antimicrobial activity against a broad spectrum of microbial organisms for up to 7 days¹
- Reduces bacterial colonization associated with catheter related infections
- Absorbs moisture around the catheter insertion site
- Silver layer is wound or skin facing with a non-adherent foam backing



| Product # | Description | Packaging | Quantity |
|-----------|-------------------------------------|-------------------------|----------|
| 46-IV20 | ¾" Disc w/ 2 mm Opening | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |
| 46-IV22 | 1" Disc w/ 4 mm Opening | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |
| 46-IV25 | 1" Disc w/ 7 mm Opening | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |
| 46-IV32 | 1" Disc w/ 4 mm Opening (w/ Insert) | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |
| 46-IV375 | 1½" Disc w/ 7 mm Opening | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |
| 46-IV34 | ¾" Disc w/out Opening | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs |


¹ Broad spectrum represented in testing by the following microbial organisms: Escherichia coli, Pseudomonas aeruginosa, Staphylococcus aureus, methicillin-resistant Staphylococcus aureus (MRSA), and Candida albicans

Covaderm® Ag+ Composite Dressings

- Provides quick and sustained antimicrobial activity against a broad spectrum of microbial organisms for up to 7 days¹
- Helps to maintain moist environment with protection from outside contaminants
- Multi-layer construction includes:
 - Protective non-adherent antimicrobial contact layer
 - Soft pad for drainage absorption
 - Conformable adhesive-coated fabric tape
- Easy application and removal
- Ideal for many wound types



Covaderm® Ag+ Composite Dressings



Covaderm® Ag+ Vascular Access Dressings (VAD)

Covaderm® Ag+ Bordered Composite Dressings

| Product # | Overall | Pad | Packaging | Quantity |
|------------|----------|-----------|-------------------------|----------|
| 46-0044AG* | 4" x 4" | 2½" x 2½" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs |
| 46-0046AG* | 4" x 6" | 2½" x 4" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs |
| 46-0066AG* | 6" x 6" | 4" x 4" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |
| 46-0048AG* | 4" x 8" | 2½" x 5½" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |
| 46-0410AG* | 4" x 10" | 2" x 7½" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |
| 46-0414AG* | 4" x 14" | 2" x 11" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |

Covaderm® Ag+ Vascular Access Dressings (VAD)

| Product # | Overall | Pad | V-Tape | Packaging | Quantity |
|-------------|---------|-----------|---------|-------------------------|----------|
| 46-0044AGV* | 4" x 4" | 2½" x 2½" | 1" x 4" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs |
| 46-0066AGV* | 6" x 6" | 4" x 4" | 1" x 6" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs |

¹ Broad spectrum represented in testing by the following microbial organisms: Escherichia coli, Pseudomonas aeruginosa, Staphylococcus aureus, methicillin-resistant Staphylococcus aureus (MRSA), and Candida albicans



***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Specialty Wound Technology Dressings

Multidex®

Multidex®

Multidex® is a hydrophilic maltodextrin wound dressing clinically proven to maintain a moist environment beneficial for granulation tissue growth and epithelial proliferation.² Multidex® fills the wound site, rapidly mixing with the exudate to form a protective coating that maintains an ideal moisture balance to protect against dehydration.

Multidex® Maltodextrin Wound Dressings

Multidex® dressing is a sterile maltodextrin composition with 1% ascorbic acid that supports an optimal wound healing environment across a variety of wound types. Available in gel or powder options, Multidex® Gel dressing is ideal for use on low to moderately exuding wounds. Multidex® Powder dressing is ideal for use in highly exuding wounds.

- Rapidly mixes with exudate within the wound to form a protective coating
- Maintains an ideal moisture balance to protect against dehydration
- For use on most wound types, including infected or hard-to-heal wounds
- Non-toxic and not systemically absorbed
- Helps control odor while decreasing exudate
- Penetrates all wound irregularities to fill tunneling and undermining
- Softens necrotic tissue to facilitate debridement

Multidex® Maltodextrin Gel Wound Dressings

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|----------------|----------|-------|
| 46-711 | ¼ fl oz | 10/Bx, 3 Bx/Cs | 30/Cs | A6248 |
| 46-710 | ½ fl oz | 10/Bx, 3 Bx/Cs | 30/Cs | A6248 |
| 46-713 | 1 fl oz | 3/Bx, 10 Bx/Cs | 30/Cs | A6248 |
| 46-712 | 3 fl oz | 6/Bx, 6 Bx/Cs | 36/Cs | A6248 |



Multidex® Maltodextrin Powder Wound Dressings

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|---------------|----------|-------|
| 46-703 | 12 g | 5/Bx, 6 Bx/Cs | 30/Cs | A6262 |
| 46-702 | 25 g | 5/Bx, 6 Bx/Cs | 30/Cs | A6262 |
| 46-701 | 45 g | 6/Bx, 6 Bx/Cs | 36/Cs | A6262 |



Foams – Absorbent Multi-Layer

Polyderm™

Polyderm™ Silicone Foam Dressings



Polyderm™ dressings are non-adherent to wounds, highly absorbent, polyurethane foam dressings.

- Ideal for moderate to highly exudative wounds
- Thick, soft construction is conformable and comfortable
- Suitable for wounds requiring protection and pressure redistribution
- Available in a variety of shapes, sizes and with and without adhesive borders

Foams – Absorbent Multi-Layer

Polyderm™

The Pressure Injury Problem:



2.5 Million people affected annually. Elderly are the most vulnerable.³



Hospital-acquired pressure injuries cost an estimated \$10,708 per patient, resulting in **\$26.8 billion spent in the US each year.**³



Estimated 60,000 individuals die annually as a direct result of pressure injuries higher mortality rate than car accidents, which kill 42,000 annually.⁴



No statistic can capture the reduction of **quality of life** such as severe pain, loss of independence, embarrassment, and increased hospital stay.

**Exceptional Care,
Accessible Price**

Polyderm™ GTL⁵ Silicone Foam Dressings

Breathable Protective Layer

The polyurethane adhesive layer glides across external surfaces, protecting wound from shear damage while allowing for moisture to escape

Superabsorbent Layer

Cushioned hydrophilic layer draws and locks exudate away from wound surface and protects surrounding periwound from maceration

Spreading Layer


Nonwoven fabric assists in distributing fluid throughout the full area of the dressing to achieve maximum absorptive capacity

Foam Layer

High-density hydrophilic foam wicks exudate away from wound surface and helps to disperse localized pressure across entire wound dressing

Perforated Silicone Layer

Gentle and atraumatic, this contact layer conforms to the anatomical contours of the skin, supplying elastic protection from shear force



Foams – Absorbent Multi-Layer

Polyderm™

Polyderm™ Silicone Foam Dressings

Polyderm™ GTL⁵ Bordered Dressings

| Product # | Description | Packaging | Quantity |
|-----------|---------------------------------------|----------------|----------|
| 46-922* | 5-Layer Bordered 2" x 2" / 5 x 5 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-933* | 5-Layer Bordered 3" x 3" / 7½ x 7½ cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-944* | 5-Layer Bordered 4" x 4" / 10 x 10 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-966* | 5-Layer Bordered 6" x 6" / 15 x 15 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



Polyderm™ GTLite Non-Bordered Dressings

| Product # | Description | Packaging | Quantity |
|-----------|---|----------------|----------|
| 46-944N* | 4-Layer Non-Bordered 4" x 4" / 10 x 10 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-966N* | 4-Layer Non-Bordered 6" x 6" / 10 x 10 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



Polyderm™ GTL⁵ Sacral Dressings

| Product # | Description | Packaging | Quantity |
|-----------|--|----------------|----------|
| 46-977* | 5-Layer Bordered Sacral Small 7" x 7" / 18 x 17 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-999* | 5-Layer Bordered Sacral Large 9" x 9" / 23 x 22 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



Polyderm™ GTL⁵ Heel Dressings

| Product # | Description | Packaging | Quantity |
|-----------|--|----------------|----------|
| 46-979* | 5-Layer Bordered Heel Large 7" x 9" / 18 x 24 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



***WARNING -** This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

DeRoyal Industries, Inc. • 200 DeBusk Lane • Powell, TN 37849 • 800.251.9864 • DeRoyal.com

Foams – Absorbent Multi-Layer

Polyderm™

Polyderm™ GTL⁵ Oval Bordered Dressings

| Product # | Description | Packaging | Quantity |
|-----------|---|----------------|----------|
| 46-924* | 5-Layer Bordered Oval 2" x 4" / 5 x 10 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



Polyderm™ GTL⁵ Dressings


| Product # | Description | Packaging | Quantity |
|-----------|--|----------------|----------|
| 46-9410* | 5-Layer Bordered 4" x 10" / 10 x 25 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



Polyderm™ GTLite Bordered Dressings

| Product # | Description | Packaging | Quantity |
|-----------|---------------------------------------|----------------|----------|
| 46-922L* | 4-Layer Bordered 2" x 2" / 5 x 5 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-944L* | 4-Layer Bordered 4" x 4" / 10 x 10 cm | 1/Pk, 10 Pk/Bx | 10/Bx |
| 46-966L* | 4-Layer Bordered 6" x 6" / 15 x 15 cm | 1/Pk, 10 Pk/Bx | 10/Bx |



 ***WARNING -** This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Foams – Absorbent Multi-Layer

Polyderm™

Polyderm™ Hydrophilic Foam Dressings

- Non-adherent, highly absorbent, lint-free polyurethane foam
- Thick, soft construction that is conformable and comfortable
- Protection from external contaminants
- Available in various styles and shapes



Polyderm™ Non-Bordered Dressings

| Product # | Overall | Foam | Packaging | Quantity | HCPCS |
|-----------|-----------|-----------|-------------------------|----------|-------|
| 46-906 | 3¾" x 3¾" | 3¾" x 3¾" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6209 |
| 46-907 | 6" x 6" | 6" x 6" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6209 |

Polyderm™ Bordered / Radial Slits Dressings

| Product # | Overall | Foam | Packaging | Quantity | HCPCS |
|-----------|-------------|-----------------------------------|-------------------------|----------|-------|
| 46-915* | 4" x 4" | 2¼" x 2¼" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6212 |
| 46-916* | 6" x 6" | 3¾" x 3¾" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6212 |
| 46-908* | 4" Diameter | 2½" Diameter w/ 2" Radial Slit | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6212 |

Polyderm™ Plus High Absorbency Foam Wound Dressings

- Polyderm™ Plus provides additional benefits to Polyderm™ Bordered dressings with the incorporation of a breathable barrier film
- Higher moisture vapor transmission rate



Polyderm™ Plus Bordered Dressings

| Product # | Overall | Foam | Packaging | Quantity | HCPCS |
|-----------|---------|-----------|-------------------------|----------|-------|
| 46-925* | 4" x 4" | 2¼" x 2¼" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6212 |
| 46-926* | 6" x 6" | 3¾" x 3¾" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6212 |



***WARNING -** This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Composite Island Dressing

Covaderm® Composite Dressings

Covaderm® Composite Dressings

Covaderm® composite dressings combine a non-adherent absorptive pad bordered by adhesive that can be applied to a variety of wound types and come in a variety of sizes and absorbency strength.

Covaderm® Composite Dressings

The Covaderm® family of general purpose sterile dressings are available in a variety of shapes, sizes, and configurations.



Composite Island Dressing

Covaderm® Composite Dressings

Covaderm® Adhesive Wound Dressings

- General purpose protective sterile dressing
- Non-adherent, absorbent sterile pad
- Non-irritating acrylic adhesive tape
- Comfortable rounded edges and stretchable fabric tape
- Split release liner for easy 1-step application
- Conforms securely to body contours



| Product # | Overall | Pad | Packaging | Quantity | HCPCS |
|-----------|----------|-----------|----------------------|----------|-------|
| 46-001* | 4" x 4" | 2½" x 2½" | 1/Pk, 25/Bx, 8 Bx/Cs | 200/Cs | A6254 |
| 46-002* | 4" x 6" | 2½" x 4" | 1/Pk, 25/Bx, 8 Bx/Cs | 200/Cs | A6254 |
| 46-003* | 4" x 8" | 2" x 5½" | 1/Pk, 25/Bx, 4 Bx/Cs | 100/Cs | A6254 |
| 46-004* | 4" x 10" | 2" x 7½" | 1/Pk, 25/Bx, 4 Bx/Cs | 100/Cs | A6254 |
| 46-005* | 4" x 14" | 2" x 11" | 1/Pk, 25/Bx, 4 Bx/Cs | 100/Cs | A6255 |

Covaderm® Thin Adhesive Wound Dressings

- Non-adherent, general purpose absorbent sterile dressing
- Conforms easily to wound
- Comfortable rounded edges and stretchable fabric tape border
- Split release liner for easy 1-step application



| Product # | Overall | Pad | Packaging | Quantity |
|-----------|----------|-----------|-----------|----------|
| 15-204 | 4" x 10" | 2½" x 7½" | 50/Cs | 50/Cs |
| 15-205 | 4" x 14" | 2½" x 11" | 50/Cs | 50/Cs |

 ***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Composite Island Dressing


Covaderm® Composite Dressings

Covaderm Plus® Adhesive Composite Dressings

- Multi-layered composite wound dressing
- Helps maintain moist environment with protection from outside contaminants
- 4-layer construction includes:
 - Protective non-adherent contact layer
 - Soft pad for drainage absorption
 - Breathable, semi-occlusive barrier film
 - Conformable adhesive-coated fabric tape
- Helps prevent drainage strike-through
- Helps minimize wound trauma
- Available with radial slit for tube sites



Covaderm Plus®
Dressings



Covaderm Plus®
Vascular Access Dressings
(VAD)

Covaderm Plus® Bordered Composite Dressings

| Product # | Overall | Pad | Packaging | Quantity | HCPCS |
|-----------|----------------------------------|--------------|-------------------------|----------|-------|
| 46-400* | 2" x 2" | 1" x 1" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-411* | 2½" x 2½" | 1½" x 1½" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6203 |
| 46-401* | 4" x 4" | 2½" x 2½" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6203 |
| 46-402* | 6" x 6" | 4" x 4" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6203 |
| 46-403* | 6" x 8" | 4" x 6" | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 50/Cs | A6204 |
| 46-414* | 4" x 10" | 2" x 7½" | 1/Pk, 100 Pk/Cs | 100/Cs | A6203 |
| 46-415* | 4" x 14" | 2" x 11" | 1/Pk, 100 Pk/Cs | 100/Cs | A6204 |
| 46-408* | 4" Diameter w/ 2" Radial Slit | 2½" Diameter | 1/Pk, 10 Pk/Bx, 5 Bx/Cs | 200/Cs | A6203 |

Covaderm Plus® Vascular Access Dressings (VAD)

| Product # | Overall | Pad | V-Tape | Packaging | Quantity | HCPCS |
|-----------|-------------|---------|---------|-------------------------|----------|-------|
| 46-405* | 4" x 4" VAD | 2" x 2" | 1" x 4" | 1/Pk, 25 Pk/Bx, 8 Bx/Cs | 200/Cs | A6254 |
| 46-406* | 6" x 6" VAD | 4" x 4" | 1" x 6" | 1/Pk, 50 Pk/Cs | 50/Cs | A6255 |



***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Contact Layers

Contact layers are non-adherent, permeable liners that conform to the wound bed, prevent trauma to fragile tissue, and allow the passage of exudate to a secondary dressing. Contact layers can also be used in combination with Negative Pressure Wound Therapy.

Dermanet® Non-Adherent Contact Layer Dressings

Dermanet® wound contact layers provide wound protection. Dermanet® and Dermanet® Ag+ (containing ionic silver) dressings are made from a proprietary blend of high-density polyethylene.



Dermanet®
Wound Contact Layers



Dermanet® Ag+ Plus
Antimicrobial Wound Contact
Layer and Absorbent Dressings



Dermanet® Ag+ / Ag+ Bordered
Antimicrobial Activity Wound
Contact Layers



Dermanet® GTL
Silicone Wound Contact Layers

Contact Layer

Dermanet® Wound Contact Layer / Dressings

Dermanet® Wound Contact Layers

- Multi-layered composite wound dressing
- Helps maintain moist environment with protection from outside contaminants
- 4-layer construction includes:
 - Protective non-adherent contact layer
 - Soft pad for drainage absorption
 - Breathable, semi-occlusive barrier film
 - Conformable adhesive-coated fabric tape
- Helps prevent drainage strike-through
- Helps minimize wound trauma



| Product # | Description | Packaging | Quantity | HCPs |
|-----------|---------------|-------------------------|----------|-------|
| 46-145 | 3" x 3" | 1/Pk, 25 Pk/Bx, 2 Bx/Cs | 50/Cs | A6206 |
| 46-143 | 5" x 4" | 1/Pk, 25 Pk/Bx, 2 Bx/Cs | 50/Cs | A6207 |
| 46-142 | 8" x 10" | 1/Pk, 25 Pk/Bx, 2 Bx/Cs | 50/Cs | A6208 |
| 46-147 | 12" x 12" | 1/Pk, 50 Pk/Cs | 50/Cs | A6208 |
| 46-141 | 20" x 15" | 1/Pk, 25 Pk/Bx, 2 Bx/Cs | 50/Cs | A6208 |
| 46-140 | 24" x 36" | 1/Pk, 50 Pk/Cs | 50/Cs | A6208 |
| 46-144 | 6" x 72" Roll | 1/Pk, 25 Pk/Cs | 25/Cs | A6208 |
| 46-G810 | 8" x 10" | 1/Pk, 25 Pk/Bx, 2 Bx/Cs | 50/Cs | A6208 |

Dermanet® Ag+ Overview

Dermanet® Ag+ dressings are a combination of Dermanet® wound contact layer and Algidex Ag®+ technology.

- A unique combination of silver alginate, maltodextrin, and Dermanet®
- Provides quick and sustained antimicrobial activity against a broad spectrum of microbial organisms for up to 7 days¹
- Does not require activation or the use of sterile distilled water
- Reduces trauma during wound dressing removal
- Decreases odor and creates a moist environment conducive to healing

Dermanet® Ag+ Contact Layers

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-DN22 | 2" x 2" | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6206 |
| 46-DN44 | 4" x 4" | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6206 |
| 46-DN48 | 4" x 8" | 1/Pk, 5 Pk/Bx, 10 Bx/Cs | 50/Cs | A6207 |
| 46-DN66 | 6" x 6" | 1/Pk, 5 Pk/Bx, 5 Bx/Cs | 25/Cs | A6207 |
| 46-DN88 | 8" x 8" | 1/Pk, 5 Pk/Bx, 5 Bx/Cs | 25/Cs | A6208 |
| 46-DN816 | 8" x 16" | 1/Pk, 5 Pk/Bx, 5 Bx/Cs | 25/Cs | A6208 |



Dermanet® Ag+ Bordered Contact Layer Dressings

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-DB22* | 2½" x 2¾" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | - |
| 46-DB44* | 4" x 4" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6206 |
| 46-DB46* | 4" x 6" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6206 |
| 46-DB410* | 4" x 10" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6207 |
| 46-DB414* | 4" x 13" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | - |



Dermanet® Ag+ Bordered Absorbent Contact Layer Dressings

| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-------------------------|----------|-------|
| 46-A44* | 4" x 4¾" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-A46* | 4" x 6" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-A410* | 4" x 10" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-A414* | 4" x 13" | 1/Pk, 25 Pk/Bx, 4 Bx/Cs | 100/Cs | A6204 |



¹ Broad spectrum represented in testing by the following microbial organisms: Escherichia coli, Pseudomonas aeruginosa, Staphylococcus aureus, methicillin-resistant Staphylococcus aureus (MRSA), and Candida albicans

WARNING - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Film Dressings

Transparent Film Dressings

Film Dressings

Transparent film dressings are adhesive, semipermeable dressings that are waterproof and impermeable to bacteria and contaminants. Films incorporate a high moisture vapor transmission rate (MVTR) design, allowing vapor and gas exchange at the wound bed.

Transeal® transparent film dressings are polyurethane films coated with an acrylic pressure-sensitive adhesive. Transeal® is highly permeable to gases and moisture vapor, allowing the skin to breathe and moisture to transmit from the site while being impermeable to water and bacteria.


Transeal®

Transparent Film Wound Dressings



Transeal® Plus

Transparent Film Wound Dressings with Absorbent Pad



Transeal® TR

Triple Release Transparent Film Wound Dressings



Transeal® MAXX

Triple Release Transparent Film Wound Dressings
(with stronger adhesive)



Transeal® Transparent Film Wound Dressings


- Transparent semi-occlusive, breathable wound dressing
- Maintains a proven high moisture vapor transmission rate⁵
- Impermeable to water and bacteria⁶
- Pressure-sensitive adhesive gently acts as second skin
- Easy, 2-step application with see-through carrier



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-----------------|----------|-------|
| 46-302* | 2½" x 2¾" | 100/Bx, 4 Bx/Cs | 400/Cs | A6257 |
| 46-303* | 4" x 4¾" | 50/Bx, 4 Bx/Cs | 200/Cs | A6258 |
| 46-306* | 8" x 12" | 10/Bx, 8 Bx/Cs | 80/Cs | A6259 |

Transeal® Plus Transparent Film Wound Dressings with Absorbent Pad

- Transparent wound dressing with non-adherent absorbent pad
- Maintains a proven high moisture vapor transmission rate⁶
- Prevents drainage strike-through
- Barrier to external contaminants



| Product # | Overall | Pad | Packaging | Quantity | HCPCS |
|-----------|----------|----------|----------------|----------|-------|
| 46-371* | 4" x 4" | 2" x 2½" | 25/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-372* | 4" x 6" | 2" x 4" | 25/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-374* | 4" x 10" | 2" x 7½" | 25/Bx, 4 Bx/Cs | 100/Cs | A6203 |
| 46-375* | 4" x 14" | 2" x 11" | 25/Bx, 4 Bx/Cs | 100/Cs | A6204 |

 ***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Film Dressings

Transparent Film Dressings

Transeal® TR Triple Release Transparent Film Wound Dressings

- Easy, three-step delivery with see-through carrier
- Transparent semi-occlusive, breathable wound dressing
- Maintains a proven high moisture vapor transmission rate⁶
- Impermeable to water and bacteria⁶
- Pressure-sensitive adhesive



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-----------|----------|-------|
| 46-325B* | 6" x 8" | 100/Cs | 100/Cs | A6258 |
| 46-326B* | 8" x 12" | 100/Cs | 100/Cs | A6258 |
| 46-302TR* | 2½" x 2¾" | 100/Cs | 100/Cs | - |

Transeal® MAXX Triple Release Transparent Film Wound Dressings


- Easy, three-step delivery with see-through carrier
- Extra strength adhesive
- Transparent semi-occlusive, breathable wound dressing
- Maintains a proven high moisture vapor transmission rate⁶
- Impermeable to water and bacteria⁶
- Pressure-sensitive adhesive



| Product # | Description | Packaging | Quantity | HCPCS |
|-----------|-------------|-----------|----------|-------|
| 46-325X | 6" x 8" | 100/Cs | 100/Cs | A6258 |
| 46-326X | 8" x 12" | 100/Cs | 100/Cs | A6259 |
| 46-302M | ½" x ¾" | 100/Cs | 100/Cs | A6257 |



***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov



The image shows a medical storage cabinet with multiple shelves. The top two shelves each contain two white DeRoyal PRO-II negative pressure wound therapy pumps. Each pump has a clear reservoir and a control panel with a digital display and buttons. Below the pumps, there are several cardboard boxes stacked on the shelves. The text "Negative Pressure Wound Therapy" is overlaid in the center of the image.

Negative Pressure Wound Therapy


Negative Pressure Wound Therapy

Continuum® NPWT Inventory Management

Continuum® NPWT System

DeRoyal's Continuum® Negative Pressure Wound Therapy (NPWT) solution aims to provide the highest level of patient care and customer satisfaction through the application of proprietary advanced clinical technologies. Combining RFID and GPS automation, DeRoyal's paperless Continuum® system allows for real-time tracking of NPWT devices, EMR integration, and inventory management that results in a quicker and safer patient discharge to home.⁷ These technologies combined with DeRoyal's powerful, simple to use, low profile negative pressure therapy unit create optimal, efficient patient outcomes.

CONTINUUM®
NPWT



One Patient, One Pump.

The DeRoyal® Negative Pressure Wound Therapy (NPWT) system, enabled with Continuum® Radio Frequency Identification (RFID) and GPS technology, is a state-of-the-art real-time supply chain management system providing anytime / anywhere dashboard access to identify negative pressure device location.

Use the NPWT Continuum® system to streamline the discharge process and enhance supply chain logistics.

Negative Pressure Wound Therapy

Continuum® NPWT Inventory Management

SCIENTIFIC/CLINICAL Benefits of VPT® technology

Variable Pressure Therapy (VPT®) technology is proven to be the most effective form of NPWT, providing a simple, comfortable and efficient wound treatment system.⁸

CARE + REVENUE Improving Care While Reducing Cost

The DeRoyal® NPWT system and protocols, enabled with the Continuum® inventory management system, ensure optimal outcomes at reduced cost levels.


LOCATE & TRACK GPS & RFID Technology

State of the art tracking systems provides detailed patient reporting, real-time supply chain management, and anytime/anywhere dashboard access to identify negative pressure device location.

INCREASED Continuity of Care

Provide continuity of care with DeRoyal's breadth and depth of product line offerings, software and hardware innovations, and devices that can seamlessly transition from hospital to home.

Hospital to DME Process



Negative Pressure Wound Therapy

Traditional Prospera® NPWT

Prospera® PRO-II™ NPWT System

The Prospera® PRO-II™ system combines powerful negative pressure therapy with hassle-free portability and clinical ease of use. Near silent operation, large, easy-to-read, color-coded display screens, and a long battery life helps increase comfort and satisfaction for patients. Offering both VPT® (Variable Pressure Technology) and Continuous settings, the Prospera® PRO-II™ unit is designed to provide personalized therapy treatment that can be used for a wide variety of wound etiologies.

Prospera® PRO-II™ Negative Pressure Wound Therapy Device

- Promotes wound healing by removal of wound exudate, debris, and infectious materials^{9,10}
- Promotes more effective blood flow to wound site^{9,10,11}
- Controls the zone of hypo-perfusion surrounding the wound edge eliminating concern for ischemia¹⁰
- Increases proliferation of new granulation tissue^{10,11,12}
- Gradual pressure cycling (in 5 mmHg increments)⁸
- Lower than traditional pressures provides gentle tissue massage resulting in less pain and better outcomes^{10,12}




| Product # | Description | Quantity |
|-------------|---------------------|----------|
| NP-2000 | PRO-II® Unit | 1/Cs |
| NP-2000GPS* | PRO-II® Unit w/ GPS | 1/Cs |

HCPCS

E2402

E2402

Variable Pressure Therapy (VPT®)¹³



Benefits of VPT® Technology

Intermittent negative pressure increases blood to the wound site more effectively than continuous negative pressure, and proliferation of new granulation tissue is significantly higher.^{9,10} Prospera® NPWT devices with VPT® technology deliver comfortable, variable pressure therapy.

- Customizable time and pressure settings
- Calibrated to increase and decrease in comfortable increments
- Limits painful high and low pressure spikes^{10,11,12}



PRO-II™ 250cc Canister PRO-II™ 450cc Canister

Prospera® PRO-II™ Canisters

NPWT Canisters are available with and without solidifier

| Product # | Size | Solidifier | Quantity |
|-----------|--------|------------|----------|
| NP-0250* | 250 cc | YES | 10/Cs |
| NP-0250N* | 250 cc | NO | 10/Cs |
| NP-0450* | 450 cc | YES | 10/Cs |
| NP-0450N* | 450 cc | NO | 10/Cs |



***WARNING** - This product can expose you to chemicals including Methylene Chloride, which is known to the State or California to cause cancer. For more information go to www.P65Warnings.ca.gov. This product can expose you to chemicals including Ethylbenzene, which is known to the State or California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



Negative Pressure Wound Therapy

Traditional NPWT Kits & Accessories

Prospera® PRO NPWT Dressing Kits

Kits: DeRoyal® NPWT kits come in a wide variety of sizes and options with each containing a filler, dome connector, and film dressing. DeRoyal's complete dressing kits also each include a measuring tool.

Black Foam: DeRoyal's black foam is a hydrophobic open-pore reticulated foam designed to promote wound healing by removing excess exudate and infection debris, and supports wound contracture when used in combination with NPWT therapy units.

White Foam: DeRoyal's white foam is an open-cell hydrophilic foam that can be used to protect delicate structures or be applied to undermining and tunneling in a wound when using NPWT therapy.

AMD Gauze: DeRoyal's antimicrobial gauze is hydrophilic and designed to absorb and retain fluid. It is effective for treatment of wounds with low to moderate drainage.




**PROPEL®
TECHNOLOGY**

- Supports efficient fluid management and accurate negative pressure delivery across all therapeutic pressure ranges
- Technology compatible with existing DeRoyal® negative pressure devices
- Technology performance supports wound dressing seal variability
- Hydrophobic Filter filters external air, limits particulate contamination into the dressing site, and protects from external fluid (e.g. showering)



**Transeal®
MAXX or TR Drape**



Black Foam



**White
Foam**



**Antimicrobial
(AMD) Gauze**

**PROPEL®
TECHNOLOGY**



**Prospera Propel®
P Dome**

Negative Pressure Wound Therapy

Traditional NPWT Kits

Prospera® PRO NPWT Kits



Black Foam Kits w/ Transeal® TR Drape

All kits contain 1 Prospera Propeller® P Dome and 1 Ruler in addition to the items listed below.

| Kit # | Contents | Quantity |
|-----------------------|---|----------|
| NP-6001P* (Small) | 1 Transeal® Triple Release Drape 8" x 12" 1 Black Foam Square (10 cm x 8 cm x 3 cm) | 10/Cs |
| NP-6002P* (Medium) | 2 Transeal® Triple Release Drapes 8" x 12" 1 Black Foam Square (20 cm x 12½ cm x 3 cm) | 10/Cs |
| NP-6003P* (Large) | 3 Transeal® Triple Release Drapes 8" x 12" 1 Black Foam Square (25 cm x 15 cm x 3 cm) | 10/Cs |



Black Foam Kits w/ Transeal® MAXX Drape

All kits contain 1 Prospera Propeller® P Dome and 1 Ruler in addition to the items listed below.

| Kit # | Contents | Quantity |
|-----------------------|--|----------|
| NP-6006P (Small) | 1 Transeal® MAXX Drape 8" x 12" 1 Black Foam Square (10 cm x 8 cm x 3 cm) | 10/Cs |
| NP-6007P (Medium) | 2 Transeal® MAXX Drapes 8" x 12" 1 Black Foam Square (20 cm x 12.5 cm x 3 cm) | 10/Cs |
| NP-6008P (Large) | 3 Transeal® MAXX Drapes 8" x 12" 1 Black Foam Square (25 cm x 15 cm x 3 cm) | 10/Cs |
| NP-0080P (X-large) | 4 Transeal® MAXX Drapes 8" x 12" 1 Black Foam Square (45 cm x 30 cm x 3 cm) | 5/Cs |



Thin Pre-cut Black Foam Kits w/ Transeal® TR Drape

All kits contain 1 Prospera Propeller® P Dome and 1 Ruler in addition to the items listed below.

| Kit # | Contents | Quantity |
|-----------------------|--|----------|
| NP-6011P* (Small) | 1 Transeal® Triple Release Drape 8" x 12" 1 Thin Helix Cut Black Foam (10 cm x 8 cm x 1.5 cm) | 10/Cs |
| NP-6012P* (Medium) | 2 Transeal® Triple Release Drapes 8" x 12" 1 Thin Helix Cut Black Foam (20 cm x 12.5 cm x 1.5 cm) | 10/Cs |
| NP-6013P* (Large) | 3 Transeal® Triple Release Drapes 8" x 12" 1 Thin Helix Cut Black Foam (25 cm x 15 cm x 1.5 cm) | 10/Cs |

⚠️ *WARNING - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Negative Pressure Wound Therapy

Traditional NPWT Kits

Prospera® PRO NPWT Kits / Specialty Dressings

Thin Pre-cut Black Foam Kits w/ Transeal® MAXX Drape

All kits contain 1 Prospera Propel® P Dome and 1 Ruler in addition to the items listed below.


| Kit # | Contents | Quantity |
|----------------------|--|----------|
| NP-6014P (Small) | 1 Transeal® MAXX Drape 8" x 12" 1 Thin Helix Cut Black Foam (10 cm x 8 cm x 1.5 cm) | 10/Cs |
| NP-6015P (Medium) | 2 Transeal® MAXX Drapes 8" x 12" 1 Thin Helix Cut Black Foam (20 cm x 12.5 cm x 1.5 cm) | 10/Cs |
| NP-6016P (Large) | 3 Transeal® MAXX Drapes 8" x 12" 1 Thin Helix Cut Black Foam (25 cm x 15 cm x 1.5 cm) | 10/Cs |



Antimicrobial Gauze Kits w/ Transeal® TR Drape

All kits contain 1 Prospera Propel® P Dome and 1 Ruler in addition to the items listed below.

| Kit # | Contents | Quantity |
|-----------------------|--|----------|
| NP-0013P* (Small) | 1 Transeal® Triple Release Drape 8" x 12" 1 AMD Gauze Sponge (2/Pk) | 10/Cs |
| NP-0014P* (Medium) | 2 Transeal® Triple Release Drapes 8" x 12" 2 AMD Gauze Sponges (2/Pk) | 10/Cs |
| NP-0015P* (Large) | 3 Transeal® Triple Release Drapes 8" x 12" 1 Sponge Kerlix® AMD Roll | 10/Cs |



Bridge Kit

All kits contain 1 Prospera Propel® P Dome in addition to the items listed below.

| Product # | Contents | Quantity |
|-----------|--|----------|
| NP-300P* | 2 Black Foam Strips 3 Transeal® Triple Release Drapes 8" x 12" 1 Thin Helix Cut Black Foam (20 cm x 12.5 cm x 1.5 cm) 1 Y-Connector | 5/Cs |



***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Negative Pressure Wound Therapy

Traditional NPWT Components / Accessories

Prospera® Dressing Components / Accessories

Prospera Propel® Dome Connector

| Product # | Description | Quantity |
|-----------|-------------|----------|
| NP-190P | Standard | 10/Cs |



Transeal® TR / MAXX Film Dressings

| Product # | Size | Quantity |
|-----------|-----------|----------|
| 46-326B* | 8" x 12" | 100/Cs |
| 46-302TR* | 2½" x 2¾" | 100/Cs |
| 46-326X | 8" x 12" | 100/Cs |
| 46-302M | ½" x ¾" | 100/Cs |



Universal Y-Connector

| Product # | Description | Quantity |
|-----------|-----------------------|----------|
| NP-120-10 | Universal Y-Connector | 10/Cs |
| NP-120 | Universal Y-Connector | 100/Cs |



 ***WARNING** - This product can expose you to chemicals including Ethylbenzene, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Negative Pressure Wound Therapy

Traditional NPWT Components / Accessories

Prospera® Dressing Components / Accessories

White Foam

| Product # | Size | Quantity |
|-----------|-----------------|----------|
| NP-906 | Small (3" x 4") | 10/Cs |
| NP-907 | Large (4" x 6") | 10/Cs |



Carrier Bags

| Product # | Description | Quantity |
|-----------|----------------------------|----------|
| NP-2001L | Prospera® PRO-II Large Bag | 1/Cs |



AC Power Adapter

| Product # | Description | Quantity |
|-----------|-----------------------------------|----------|
| NP-2002 | Prospera® PRO-II AC Power Adapter | 1/Cs |



Bed Rack / I.V. Pole Attachments

| Product # | Description | Quantity |
|-----------|-----------------------|----------|
| NP-170 | Bed Rack | 1/Cs |
| NP-7205 | I.V. Pole | 1/Cs |
| NP-170IV | Bed Rack w/ I.V. Pole | 1/Cs |



Negative Pressure Wound Therapy

Disposable Prospera® NPWT

Prospera Flex® NPWT System

A disposable, single-use NPWT system containing a device, a canister, tubing, and a wound dressing interface. It is indicated for use on a variety of wound types, including closed surgical incisions. It can be used with dressings cleared for NPWT, providing clinicians flexibility to choose the most appropriate dressing for their patients.

NPWT Benefits Include

Acts as a Barrier to External Contamination^{14,15}



Provides an environment to promote wound healing^{14,15}



Helps to hold Incision Edges Together¹⁶



Removes Fluids and Infectious Materials^{14,15}



Decreases Lateral Tensions of Sutured/Stapled Incisions^{16,17}



Reduces Edema¹⁵



Prospera Flex® Single-Use 7 Day NPWT Device

- » Single button user interface
- » Provides 7 days of NPWT
- » Delivers -80 mmHg of continuous NPWT
- » Audible and visual alarms
- » 45 cc flexible canister
- » Canister integrated tubing (4ft)
- » Discrete carrying case
- » AA battery powered



Prospera Propel® FX Dome

- » Incorporates Prospera Propel® Technology for efficient fluid removal
- » Self-sealing dome prevents fluid from entering wound interface when tubing is not connected
- » Connector at dome allows for convenient tubing detachment, supporting patient daily activities (showering, clothing, etc.)

Prospera Flex® NPWT System is Compatible with Lumitech® Incisional Wound Dressings

- Distributes pressure evenly across the dressing
- Reduces lateral tension when used over closed incisions
- Available together in a post op convenience kit for ease of ordering and application

Negative Pressure Wound Therapy

Disposable Prospera® NPWT

Prospera Flex® NPWT System

Product #
NP-8007S*

Description
Prospera Flex® NPWT Single-Use System
Device, Canister, Dome, Accessories

Quantity
5 Ea/Cs

NP-8007S-01*

Prospera Flex® NPWT Single-Use System
Device, Canister, Dome, Accessories

1 Ea



Prospera Flex® Canisters

Product #
NP-8045C*

Description
Prospera Flex® Canister – 45 cc

Quantity
10 Ea/Cs



Prospera Flex® Dome

Product #
NP-190F*

Description
Prospera Propel® FX Dome

Quantity
5 Ea/Cs



***WARNING -** This product can expose you to chemicals including Vinylchloride, which is known to the State or California to cause cancer, and 1,2-Butadiene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

This product can expose you to chemicals including Methylene Chloride, which is known to the State or California to cause cancer, For more information go to www.P65Warnings.ca.gov

Negative Pressure Wound Therapy


NPWT Compatible Dressings & Accessories

Lumitech® Incisional Dressing


The Lumitech® multi-layer incisional wound dressing is intended to manage post-surgical wounds with low to medium exudate levels. The dressing incorporates Lumina® Adhesive Technology, a patented UV switchable adhesive that supports a strong seal and long wear time yet allows for gentle dressing removal. This helps prolong dressing usage while minimizing painful dressing removal.

Innovative Features

- » Peel and place, sterile multi-layer absorbent design
- » Integrated Lumina® Adhesive Technology
- » Optimized design for application over various incision site locations
- » NPWT compatible dressing with dome site positioning flexibility
- » Available in a variety of sizes



LUMITECH®




Multi-Layer Dressing Composition

1. Integrated wound contact layer with extended edges for increased skin integrity and patient comfort¹⁸
2. Non-woven fibers support up to 20 ccs of exudate absorption at the wound site when used independently or with NPWT¹⁹
3. Polyurethane foam assists with pressure distribution when used with NPWT¹⁹
4. Polyurethane film incorporating Lumina® Adhesive Technology provides flexible, breathable, water-resistant, contaminant-blocking protection with a strong seal and gentle removal
5. Release liners and carrier layers for ease of application

Lumina® Adhesive Technology Benefits²⁰

- » Adhesive adheres as strong as acrylic dressings, supporting long wear time in unswitched state
- » Adhesives releases more gently than silicone following UV light exposure, minimizes pain and periwound trauma at dressing removal



**Compatible with the DeRoyal® Prospera Flex®
NPWT Single-Use System**

Negative Pressure Wound Therapy

NPWT Compatible Dressings & Accessories

Incision Management Convenience Kit

| Product # | Description | Quantity |
|-------------|--|----------|
| 47-813K* | Incision Management Convenience Kit (Small) x1 Prospera Flex® NPWT System x2 Lumitech® 13 cm Incisional Wound Dressings x1 Lumitech® Compatible UV Flashlight Accessories | 5 Ea/Cs |
| 47-813K-01* | Incision Management Convenience Kit (Small) x1 Prospera Flex® NPWT System x2 Lumitech® 13 cm Incisional Wound Dressings x1 Lumitech® Compatible UV Flashlight Accessories | 1 Ea |
| 47-820K* | Incision Management Convenience Kit (Medium) x1 Prospera Flex® NPWT System x2 Lumitech® 20 cm Incisional Wound Dressings x1 Lumitech® Compatible UV Flashlight Accessories | 5 Ea/Cs |
| 47-820K-01* | Incision Management Convenience Kit (Medium) x1 Prospera Flex® NPWT System x2 Lumitech® 20 cm Incisional Wound Dressings x1 Lumitech® Compatible UV Flashlight Accessories | 1 Ea |



Lumitech® Dressings

| Product # | Description | Quantity |
|-----------|-------------------------------------|----------|
| 46-LT13 | Lumitech® 13 cm Incisional Dressing | 10 Ea/Cs |
| 46-LT20 | Lumitech® 20 cm Incisional Dressing | 10 Ea/Cs |



Lumitech® Dressing Accessories

| Product # | Description | Quantity |
|--------------|-------------------------------------|-------------------|
| 46-LTUV | Lumitech® Compatible UV Flashlight | 10 Ea/Cs |
| 46-LT20S-30 | Lumitech® Securement Strips (20 cm) | 6 Ea/Pk, 5 Pk/Cs |
| 46-LT20S-120 | Lumitech® Securement Strips (20 cm) | 6 Ea/Pk, 20 Pk/Cs |



***WARNING** - This product can expose you to chemicals including Vinylchloride, which is known to the State of California to cause cancer, and Toluene, which is known to the State of California to cause birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov

Appendix

WOUND CARE HCPCS

HCPCS²¹ information is for informational purposes only. It represents no statement, promise, or guarantee by DeRoyal concerning level of reimbursement, payment, or charge. It does not and is not intended to increase or maximize reimbursement, nor does it promise or guarantee that the codes listed will be appropriate or that reimbursement will be made. Please consult your payor organization to obtain the most current information or local Medicare Part B coverage and reimbursement policies.

| HCPCS | Long Description | Exudate | Frequency of Change | DMEPOS |
|-------|--|----------------------------|---------------------|--------|
| A4450 | Tape, non-waterproof, per 18 square inches | | NSP | YES |
| A4452 | Tape, waterproof, per 18 square inches | | NSP | YES |
| A4461 | Surgical dressing holder, non-reusable, each | | NSP | YES |
| A4463 | Surgical dressing holder, reusable, each | | NSP | YES |
| A4490 | Surgical stockings above knee length, each | | NSP | NO |
| A4495 | Surgical stockings thigh length, each | | NSP | NO |
| A4500 | Surgical stockings below knee length, each | | NSP | NO |
| A4510 | Surgical stockings full length, each | | NSP | NO |
| A4649 | Surgical supply; miscellaneous | | NSP | NO |
| A6010 | Collagen based wound filler, dry form, sterile, per gram of collagen | | NSP | YES |
| A6011 | Collagen based wound filler, gel/paste, per gram of collagen | | NSP | YES |
| A6021 | Collagen dressing, sterile, size 16 sq. in. or less, each | | NSP | YES |
| A6022 | Collagen dressing, sterile, size more than 16 sq. in. but less than or equal to 48 sq. in., each | | NSP | YES |
| A6023 | Collagen dressing, sterile, size more than 48 sq. in., each | | NSP | YES |
| A6024 | Collagen dressing wound filler, sterile, per 6 inches | | NSP | YES |
| A6025 | Gel Sheet for dermal or epidermal application, (e.g., silicone, hydrogel, other), each | | NSP | NO |
| A6154 | Wound pouch, each | | 3 times /week | YES |
| A6196 | Alginate or other fiber gelling dressing, wound cover, sterile, pad size 16 sq. in. or less, each dressing | Heavy exudation | Once daily | YES |
| A6197 | Alginate or other fiber gelling dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., each dressing | Moderate-Heavy Exudation | Once daily | YES |
| A6198 | Alginate or other gelling dressing, wound cover, sterile, pad size more than 48 sq. in., each dressing | Moderate-Heavy Exudation | Once daily | NO |
| A6199 | Alginate or other fiber gelling dressing, wound filler, sterile, per 6 inches | Moderate-Heavy Exudation | 12 inches per day | YES |
| A6203 | Composite dressing, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | | 3 times /week | YES |
| A6204 | Composite dressing, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | | 3 times /week | YES |
| A6205 | Composite dressing, sterile, pad size more than 48 sq. in. with any size adhesive border, each dressing | | 3 times /week | NO |
| A6206 | Contact layer, sterile, 16 sq. in. or less, each dressing | | Once/week | NO |
| A6207 | Contact layer, sterile, more than 16 sq. in. but less than or equal to 48 sq. in., each dressing | | Once/week | YES |
| A6208 | Contact layer, sterile, 48 sq. in. or more each dressing | | Once/week | NO |
| A6209 | Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | YES |
| A6210 | Foam dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | YES |
| A6211 | Foam dressing, wound cover, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | YES |
| A6212 | Foam dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | YES |
| A6213 | Foam dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | NO |
| A6214 | Foam dressing, wound cover, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing | Moderate-Heavy Exudation | 3 times /week | YES |
| A6215 | Foam dressing, wound filler, sterile, per gram | Moderate-Heavy Exudation | Once/week | NO |
| A6216 | Gauze, non-impregnated, non-sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | | 3/day | YES |
| A6217 | Gauze, non-impregnated, non-sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | | 3/day | YES |
| A6218 | Gauze, non-impregnated, non-sterile, pad size more than 48 sq. in., without adhesive border, each dressing | | 3/day | NO |
| A6219 | Gauze, non-impregnated, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | | Once daily | YES |
| A6220 | Gauze, non-impregnated, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | | Once daily | YES |
| A6221 | Gauze, non-impregnated, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing | | Once daily | NO |
| A6222 | Gauze, impregnated with other than water, normal saline, or hydrogel, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | | Once daily | YES |
| A6223 | Gauze, impregnated with other than water, normal saline, or hydrogel, sterile, pad size more than 16 sq. in., but less than or equal to 48 sq. in., without adhesive border, each dressing | | Once daily | YES |
| A6224 | Gauze, impregnated with other than water, normal saline, or hydrogel, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | | Once daily | YES |
| A6228 | Gauze, impregnated with water or normal saline, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | | NSP | NO |
| A6229 | Gauze, impregnated, water or normal saline, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | | NSP | YES |
| A6230 | Gauze, impregnated with water or normal saline, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | | NSP | NO |
| A6231 | Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size 16 sq. in. or less, each dressing | Minimum Exudate | Once daily | YES |
| A6232 | Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size greater than 16 sq. in., but less than or equal to 48 sq. in., each dressing | Minimum Exudate | Once daily | YES |
| A6233 | Gauze, impregnated, hydrogel, for direct wound contact, sterile, pad size more than 48 sq. in., each dressing | Minimum Exudate | Once daily | YES |
| A6234 | Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | YES |
| A6235 | Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | YES |
| A6236 | Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | YES |
| A6237 | Hydrocolloid dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | YES |
| A6238 | Hydrocolloid dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | YES |
| A6239 | Hydrocolloid dressing, wound cover, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing | Light - Moderate Exudation | 3 times /week | NO |
| A6240 | Hydrocolloid dressing, wound filler, paste, sterile, per ounce | Light - Moderate Exudation | 3 times /week | YES |
| A6241 | Hydrocolloid dressing, wound filler, dry form, sterile, per gram | Light - Moderate Exudation | 3 times /week | YES |
| A6242 | Hydrogel dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | Minimum Exudate | Once daily | YES |
| A6243 | Hydrogel dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | Minimum Exudate | Once daily | YES |
| A6244 | Hydrogel dressing, wound cover, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | Minimum Exudate | Once daily | YES |
| A6245 | Hydrogel dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | Minimum Exudate | 3 times /week | YES |
| A6246 | Hydrogel dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | Minimum Exudate | 3 times /week | YES |
| A6247 | Hydrogel dressing, wound cover, sterile, pad size more than 48 sq. in., with any size adhesive border, each dressing | Minimum Exudate | 3 times /week | YES |
| A6248 | Hydrogel dressing, wound filler, gel, per fluid ounce | Minimum Exudate | 3 oz. /month | YES |
| A6250 | Skin sealants, protectants, moisturizers, ointments, any type, any size | | NSP | NO |

WOUND CARE HCPCS

HCPCS²¹ information is for informational purposes only. It represents no statement, promise, or guarantee by DeRoyal concerning level of reimbursement, payment, or charge. It does not and is not intended to increase or maximize reimbursement, nor does it promise or guarantee that the codes listed will be appropriate or that reimbursement will be made. Please consult your payor organization to obtain the most current information or local Medicare Part B coverage and reimbursement policies.

| HCPCS | Long Description | Exudate | Frequency of Change | DMEPOS |
|-------|---|--------------------------|---------------------|--------|
| A6251 | Specialty absorptive dressing, wound cover, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | Moderate-Heavy Exudation | Once daily | YES |
| A6252 | Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., without adhesive border, each dressing | Moderate-Heavy Exudation | Once daily | YES |
| A6253 | Specialty absorptive dressing, wound cover, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | Moderate-Heavy Exudation | Once daily | YES |
| A6254 | Specialty absorptive dressing, wound cover, sterile, pad size 16 sq. in. or less, with any size adhesive border, each dressing | Moderate-Heavy Exudation | Every 48 hours | YES |
| A6255 | Specialty absorptive dressing, wound cover, sterile, pad size more than 16 sq. in. but less than or equal to 48 sq. in., with any size adhesive border, each dressing | Moderate-Heavy Exudation | Every 48 hours | YES |
| A6256 | Specialty absorptive dressing, wound cover, sterile, more than 48 sq. in., with any size adhesive border, each dressing | Moderate-Heavy Exudation | Every 48 hours | NO |
| A6257 | Transparent film, sterile, 16 sq. in. or less, each dressing | Minimum Exudate | 3 times /week | YES |
| A6258 | Transparent film, sterile, more than 16 sq. in. but less than or equal to 48 sq. in., each dressing | Minimum Exudate | 3 times /week | YES |
| A6259 | Transparent film, sterile, more than 48 sq. in., each dressing | Minimum Exudate | 3 times /week | YES |
| A6260 | Wound cleanser, any type, any size | | NSP | NO |
| A6261 | Wound filler, gel/paste, per liquid ounce, not otherwise classified | | Once daily | NO |
| A6262 | Wound filler, gel/paste, per gram, not otherwise classified | | Once daily | NO |
| A6266 | Gauze, impregnated, other than water, normal saline, or zinc paste, sterile, any width, per linear yard | | Once daily | YES |
| A6402 | Gauze, non-impregnated, sterile, pad size 16 sq. in. or less, without adhesive border, each dressing | | 3 times/day | YES |
| A6403 | Gauze, non-impregnated, sterile, pad size more than 16 sq. in. less than or equal to 48 sq. in., without adhesive border, each dressing | | 3 times/day | YES |
| A6404 | Gauze, non-impregnated, sterile, pad size more than 48 sq. in., without adhesive border, each dressing | | 3 times/day | NO |
| A6407 | Packing strips, non-impregnated, sterile, up to 2 inches in width, per linear yard | | 3 times/day | YES |
| A6410 | Eye pad, sterile, each | | NSP | YES |
| A6411 | Eye pad, non-sterile, each | | NSP | YES |
| A6412 | Eye patch, occlusive, each | | NSP | NO |
| A6413 | Adhesive bandage, first-aid type, any size each | | NSP | NO |
| A6441 | Padding bandage, non-elastic, non-woven/non-knitted, width greater than or equal to three inches and less than five inches, per yard | | Frequency varies | YES |
| A6442 | Conforming bandage, non-elastic, knitted/woven, non-sterile, width less than three inches, per yard | | Frequency varies | YES |
| A6443 | Conforming bandage, non-elastic, knitted/woven, non-sterile, width greater than or equal to three inches and less than five inches, per yard | | Frequency varies | YES |
| A6444 | Conforming bandage, non-elastic, knitted/woven, non-sterile, width greater than or equal to 5 inches, per yard | | Frequency varies | YES |
| A6445 | Conforming bandage, non-elastic, knitted/woven, sterile, width less than three inches, per yard | | Frequency varies | YES |
| A6446 | Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to three inches and less than five inches, per yard | | Frequency varies | YES |
| A6447 | Conforming bandage, non-elastic, knitted/woven, sterile, width greater than or equal to five inches, per yard | | Frequency varies | YES |
| A6448 | Light compression bandage, elastic, knitted/woven, width less than three inches, per yard | | Once /week | YES |
| A6449 | Light compression bandage, elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard | | Once /week | YES |
| A6450 | Light compression bandage, elastic, knitted/woven, width greater than or equal to five inches, per yard | | Once /week | YES |
| A6451 | Moderate compression bandage, elastic, knitted/woven, load resistance of 1.25 to 1.34 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard | | Once /week | YES |
| A6452 | High compression bandage, elastic, knitted/woven, load resistance greater than or equal to 1.35 foot pounds at 50% maximum stretch, width greater than or equal to three inches and less than five inches, per yard | | Once /week | YES |
| A6453 | Self-adherent bandage, elastic, non-knitted/non-woven, width less than three inches, per yard | | Once /week | YES |
| A6454 | Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to three inches and less than five inches, per yard | | Once /week | YES |
| A6455 | Self-adherent bandage, elastic, non-knitted/non-woven, width greater than or equal to five inches, per yard | | Once /week | YES |
| A6456 | Zinc paste impregnated bandage, non-elastic, knitted/woven, width greater than or equal to three inches and less than five inches, per yard | | Once /week | YES |
| A6457 | Tubular dressing with or without elastic, any width, per linear yard | | NSP | YES |
| A6501 | Compression burn garment, bodysuit (head to foot), custom fabricated | | NSP | YES |
| A6502 | Compression burn garment, chin strap, custom fabricated | | NSP | YES |
| A6503 | Compression burn garment, facial hood, custom fabricated | | NSP | YES |
| A6504 | Compression burn garment, glove to wrist, custom fabricated | | NSP | YES |
| A6505 | Compression burn garment, glove to elbow, custom fabricated | | NSP | YES |
| A6506 | Compression burn garment, glove to axilla, custom fabricated | | NSP | YES |
| A6507 | Compression burn garment, foot to knee length, custom fabricated | | NSP | YES |
| A6508 | Compression burn garment, foot to thigh length, custom fabricated | | NSP | YES |
| A6509 | Compression burn garment, upper trunk to waist including arm openings (vest), custom fabricated | | NSP | YES |
| A6510 | Compression burn garment, trunk, including arms down to leg openings (leotard), custom fabricated | | NSP | YES |
| A6511 | Compression burn garment, lower trunk including leg openings (panty), custom fabricated | | NSP | YES |
| A6512 | Compression burn garment, not otherwise classified | | NSP | NO |
| A6513 | Compression burn mask, face and/or neck, plastic or equal, custom fabricated | | NSP | YES |
| A6530 | Gradient compression stocking, below open venous ulcer knee, 18-30 mm HG, each | | NSP | NO |
| A6531 | Gradient compression stocking, below knee, 30-40 mmhg, each | | NSP | YES |
| A6532 | Gradient compression stocking, below knee, 40-50 mmhg, each | | NSP | YES |
| A6533 | Gradient compression stocking, thigh open venous ulcer length, 18-30 mm HG, each | | NSP | NO |
| A6534 | Gradient compression stocking, thigh open venous ulcer length, 30-40 mm HG, each | | NSP | NO |
| A6535 | Gradient compression stocking, thigh open venous ulcer length, 40-50 mm HG, each | | NSP | NO |
| A6536 | Gradient compression stocking, full open venous ulcer length/chap style, 18-30 mm HG, each | | NSP | NO |
| A6537 | Gradient compression stocking, full open venous ulcer length/chap style, 30-40 mm HG, each | | NSP | NO |
| A6538 | Gradient compression stocking, full open venous ulcer length/chap style, 40-50 mm HG, each | | NSP | NO |
| A6539 | Gradient compression stocking, waist open venous ulcer length, 18-30 mm HG, each | | NSP | NO |
| A6540 | Gradient compression stocking, waist open venous ulcer length, 30-40 mm HG, each | | NSP | NO |
| A6541 | Gradient compression stocking, waist open venous ulcer length, 40-50 mm HG, each | | NSP | NO |
| A6544 | Gradient compression stocking, garter belt | | NSP | NO |
| A6545 | Gradient compression wrap, non-elastic, below knee, 30-50 mm hg, each | | NSP | YES |
| A9720 | Non-covered item or service | | NSP | NO |

NSP = No standard provided by Medicare, refer to Medicare utilization guidelines²²

Frequency Varies = Frequency and approval of medical necessity will vary based on medical condition, primary or secondary dressing choices and wound exudation documentation²¹

DMEPOS = If code is listed in 2014 DMEPOS, Medicare has pre-determined the reimbursement for this item, refer to guidelines²³

Attributions & References

Trademark Attributions

DeRoyal, the DeRoyal logo, Algidex Ag, Apex, Continuum, the Continuum Logo, Covaderm, Covaderm Plus, Dermanet, DeWrap, Episeal, Flutext, Lumitech, Multidex, Multipad, Polyderm, Prospera Flex, Prospera Propel, PRO-II, Samson, Sofisorb, StayDry, Stretch Net, Transeal, and VPT are trademarks or registered trademarks of DeRoyal Industries, Inc. in the United States and other countries.

JetOx trademark is used under license by DeRoyal Industries, Inc.

Kerlix is a registered trademark of KPR U.S., LLC

Lumina is a registered trademark of Lumina Adhesives AB.

References

1. Data on File.
2. Data on File.
3. Padula; Delarmente. The national cost of hospital-acquired pressure injuries in the United States. *Int Wound J.* 2019 Jun;16(3):634–640. doi: 10.1111/ iwj.13071. Epub 2019 Jan 28.
4. US Dept. of Transportation, National Highway Traffic Safety Administration, 2022.
5. Data on File.
6. Data on File.
7. Jensen R, Cayce J, Galinis J (2019). Improving Quality of Care and Reducing Cost Burdens Through NPWT Management Automation. Society of Advanced Wound Care, San Antonio, TX.
8. Malmjö, M., Gustafsson, L., Lindstedt, S., Gesslein, B., & Ingemansson, R. (2012). The Effects of Variable, Intermittent, and Continuous Negative Pressure Wound Therapy, Using Foam or Gauze, on Wound Contraction, Granulation Tissue Formation, and Ingrowth Into the Wound Filler. *Eplasty*, 12, e5.
9. Morykwas MJ, Argenta LC, Shelton-Brown EI, McGuirt W. Vacuum-Assisted Closure: A new method for wound control and treatment: Animal studies and basic foundation. *Annals of Plastic Surgery* 1997;38:553-562.
10. Borquist, O., Ingemansson, R., Malmjö, M., The Influence of Low and High Pressure Levels during Negative Pressure Wound Therapy on Wound Contraction and Fluid Evacuation. *Plastic and Reconstructive Surgery*, February 2011; 127:2: 551-559.
11. Eberlein T, Fendler H, Ahearn C. Using a new technique for negative pressure wound therapy (NPWT) for the management of chronic, non-healing wounds. *SAWC 2008*, San Diego.
12. Eberlein T, Fendler H. Using a new technique of negative pressure wound therapy (NPWT), Variable Pressure Therapy (VPT) for the management of chronic, non-healing wounds. *JWOCN* May/June 2009;36(35):S23.
13. Data on File.
14. Normandin et al. Negative Pressure Wound Therapy: Mechanism of Action and Clinical Applications. *Semin Plast Surg.* 2021 Aug; 35(3): 164–170.
15. Quacinella et al. Negative pressure wound therapy: Where are we in 2022? *OTA Int.* 2023 Jul; 6(4 Suppl): e247.
16. Dohmen et al. Use of Incisional Negative Pressure Wound Therapy on Closed Median Sternal Incisions after Cardiothoracic Surgery: Clinical Evidence and Consensus *Med Sci Monit.* 2014; 20: 1814–1825.
17. Data on File.
18. Stocks et al. "Evaluation of a Novel Incisional Dressing and Disposable Negative Pressure Device Compared to the Market-Leading System for Dressing Survivability, User Comfort, and Skin Cleanliness: A randomized trial in Healthy Volunteers" *SAWC Fall 2024*; 10.13140/RG.2.2.18141.83686.
19. Data on File.
20. Cayce, Jonathan, Ph.D., MS, Goins, Lily, & Warrington, Gavin. (2023, April 26). Light-Deactivated Adhesive: The Solution to Maintaining Dressing Integrity and Protecting Skin. Society for Advanced Wound Care, National Harbor, MD.
21. Centers for Medicare & Medicaid Services. List of CPT/HCPCS Codes | CMS. Retrieved May 28, 2028, from <https://www.cms.gov/medicare/regulations-guidance/physician-self-referral/list-cpt-hcpcs-codes>.
22. Centers for Medicare & Medicaid Services. Medicare Fee for Service for Parts A & B | CMS. Retrieved May 28, 2028, from <https://www.cms.gov/data-research/statistics-trends-and-reports/medicare-fee-for-service-parts-a-b>.
23. Centers for Medicare & Medicaid Services. Durable Medical Equipment, Prosthetic Devices, Prosthetics, Orthotics, & Supplies | CMS. Retrieved May 28, 2028, from <https://www.cms.gov/medicare/payment/fee-schedules/dmepos>.



DeRoyal Industries, Inc.

200 DeBusk Lane, Powell, TN 37849 USA

888.938.7828 | DeRoyal.com