

# Sizing Guide

Sizing determined by three circumferential measurements:

1. 1" below tip of breast bone or xyphoid process
2. Around navel or waist
3. Around fullest area of trochanter or hips

## Non-Tapered

Small

Medium

Large

X-Large

XX-Large

## Xyphoid

31" - 35"

34" - 38"

37" - 41"

38" - 45"

45" - 53"

## Waist

31" - 35"

32" - 38"

35" - 41"

38" - 45"

45" - 53"

## Trochanter

31" - 35"

34" - 39"

38" - 42"

38" - 45"

45" - 53"

## Tapered

Small

Medium

Large

X-Large

XX-Large

## Xyphoid

25" - 31"

28" - 35"

32" - 38"

35" - 42"

40" - 46"

## Waist

25" - 31"

29" - 36"

32" - 38"

36" - 42"

39" - 47"

## Trochanter

27" - 33"

32" - 38"

35" - 41"

38" - 46"

43" - 49"